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

11 Bibliographie

1. Y. Okada, H. Nagase, and E. D. Harris, Jr. A metalloproteinase from human rheumatoid synovial fibroblasts that digests connective tissue matrix components. Purification and characterization. *J.Biol.Chem.* 261 (30):14245-14255, 1986.
2. K. Obata, K. Iwata, Y. Okada, Y. Kohrin, E. Ohuchi, S. Yoshida, M. Shinmei, and T. Hayakawa. A one-step sandwich enzyme immunoassay for human matrix metalloproteinase 3 (stromelysin-1) using monoclonal antibodies. *Clin.Chim.Acta* 211 (1-2):59-72, 1992.
3. J. Martel-Pelletier, R. McCollum, N. Fujimoto, K. Obata, JM. Cloutier, JP. Pelletier. Excess of metalloproteinases over tissue inhibitor of metalloproteinase may contribute to cartilage degradation in osteoarthritis and rheumatoid arthritis. *Lab Invest.* 70 (6):807-15, 1994.
4. S. Sasaki, H. Iwata, N. Ishiguro, K. Obata, and T. Miura. Detection of stromelysin in synovial fluid and serum from patients with rheumatoid arthritis and osteoarthritis. *Clin.Rheumatol.* 13 (2):228-233, 1994.
5. Y. Yoshihara, K. Obata, N. Fujimoto, K. Yamashita, T. Hayakawa, and M. Shimmei. Increased levels of stromelysin-1 and tissue inhibitor of metalloproteinases-1 in sera from patients with rheumatoid arthritis. *Arthritis Rheum.* 38 (7):969-975, 1995.
6. H. Yamanaka, Y. Matsuda, M. Tanaka, W. Sendo, H. Nakajima, A. Taniguchi, and N. Kamatani. Serum matrix metalloproteinase 3 as a predictor of the degree of joint destruction during the six months after measurement, in patients with early rheumatoid arthritis. *Arthritis Rheum.* 43 (4):852-858, 2000.
7. M. D. Posthumus, P. C. Limburg, J. Westra, M. A. Van Leeuwen, and M. H. van Rijswijk. Serum matrix metalloproteinase 3 levels during treatment with sulfasalazine or combination of methotrexate and sulfasalazine in patients with early rheumatoid arthritis. *J.Rheumatol.* 29 (5):883-889, 2002.
8. Katrib, M. D. Smith, M. J. Ahern, J. Slavotinek, L. Stafford, C. Cuello, J. V. Bertouch, H. P. McNeil, and P. P. Youssef. Reduced chemokine and matrix metalloproteinase expression in patients with rheumatoid arthritis achieving remission. *J.Rheumatol.* 30 (1):10-21, 2003.
9. Tchvetverikov, L. R. Lard, J. DeGroot, N. Verzijl, J. M. TeKoppele, F. C. Breedveld, T. W. Huizinga, and R. Hanemaaijer. Matrix metalloproteinases-3, -8, -9 as markers of disease activity and joint damage progression in early rheumatoid arthritis. *Ann.Rheum.Dis.* 62 (11):1094-1099, 2003.
10. M. J. Green, A. K. Gough, J. Devlin, J. Smith, P. Astin, D. Taylor, and P. Emery. Serum MMP-3 and MMP-1 and progression of joint damage in early rheumatoid arthritis. *Rheumatology.(Oxford)* 42 (1):83-88, 2003.
11. Litinsky, D. Paran, D. Levartovsky, I. Wigler, I. Kaufman, I. Yaron, M. Yaron, D. Caspi, O. Elkayam. The effects of leflunomide on clinical parameters and serum levels of IL-6, IL-10, MMP-1 and MMP-3 in patients with resistant rheumatoid arthritis. *Cytokine.* 33 (2):106-10, 2006.
12. S. Visvanathan, M. U. Rahman, E. Keystone, M. Genovese, L. Klareskog, E. Hsia, M. Mack, J. Buchanan, M. Elashoff, and C. Wagner. Association of serum markers with improvement in clinical response measures after treatment with golimumab in patients with active rheumatoid arthritis despite receiving methotrexate: results from the FORWARD study. *Arthritis Res.Ther.* 12 (6):R211, 2010.
13. Mahemara, T. Sugimoto, D. Sugiyama, S. Morinobu, G. Tsuji, S. Kawano, A. Morinobu, and S. Kumagai. Serum Matrix Metalloproteinase-3 as Predictor of Joint Destruction in Rheumatoid Arthritis, Treated with Non-biological Disease Modifying Anti-Rheumatic Drugs. *Kobe J. Med. Sci.* 56 (3):98-107, 2010.
14. Y. Urata, R. Uesato, D. Tanaka, Y. Nakamura, and S. Motomura. Treating to target matrix metalloproteinase 3 normalisation together with disease activity score below 2.6 yields better effects than each alone in rheumatoid arthritis patients: T-4 Study. *Ann.Rheum.Dis.*, 2011.
15. N. Nishimoto, K. Amano, Y. Hirabayashi, T. Horiuchi, T. Ishii, M. Iwahashi, M. Iwamoto, H. Kohsaka, M. Kondo, T. Matsubara, T. Mimura, H. Miyahara, S. Ohta, Y. Saeki, K. Saito, H. Sano, K. Takasugi, T. Takeuchi, S. Tohma, T. Tsuru, Y. Ueki, J. Yamana, J. Hashimoto, T. Matsutani, M. Murakami, and N. Takagi. Drug free REmission/low disease activity after cessation of tocilizumab (Actemra) Monotherapy (DREAM) study. *Mod.Rheumatol.*, 2013.
16. Hiura, S. Iwaki-Egawa, T. Kawashima, S. I. Fujisawa, T. Takeda, H. Komori, and Y. Watanabe. The diagnostic utility of matrix metalloproteinase-3 and high-sensitivity C-reactive protein for predicting rheumatoid arthritis in anti-cyclic citrullinated peptide antibody-negative patients with recent-onset undifferentiated arthritis. *Rheumatol.Int.*, 2013.

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IVD	- Diagnosi in vitro - Pour diagnostic in vitro - In Vitro Diagnostikum - Para uso Diagnóstico in vitro	- For in vitro diagnostic use - Para uso diagnóstico in vitro - In Vitro Διαγνωστικό μέσο
REF	- Numero d'ordine - Référence Catalogue - Bestellnummer - Número de catálogo	- Catalogue number - Numéro de catálogo - Αριθμός παραγγελίας
LOT	- Descrizione lotto - Lot - Chargen Bezeichnung - Lote	- Lot - Lote - Χαρακτηρισμός παρτίδας
UDI	- Identificatore univoco del dispositivo - Identifiant unique de l'appareil - eindeutige Produktidentifizierung - Identificador único do dispositivo	- Unique Device Identifier - Identificador único del dispositivo - Μοναδικό αναγνωριστικό συσκευής
CE	- Conformità europea - Déclaration CE de Conformité - Europäische Konformität - Declaração CE de Conformidade	- EC Declaration of Conformity - Declaración CE de Conformidad - Ευρωπαϊκή συμφωνία
	- 96 determinazioni - 96 tests - 96 Bestimmungen - 96 Testes	- 96 tests - 96 pruebas - 96 προσδιορισμοί
	- Rispettare le istruzioni elettroniche per l'uso - Voir les instructions d'utilisation électronique - Elektronische Gebrauchsanweisung beachten - Seguir as instruções electrónicas de utilização	- See electrical instructions for use - Sigas las instrucciones electrónicas de uso - Ακολουθήστε τις ηλεκτρονικές οδηγίες χρήσης
	- Da utilizzarsi entro - Utilise avant le - Verwendbar bis - Utilizar antes de	- Use by - Utilizar antes de - Χρήση μέχρι
	- Conservare a 2-8°C (35,6-46,4°F) - Conserver à 2-8°C (35,6-46,4°F) - Lagerung bei 2-8°C (35,6-46,4°F) - Conservar entre 2-8°C (35,6-46,4°F)	- Store at 2-8°C (35,6-46,4°F) - Conservar a 2-8°C (35,6-46,4°F) - Φυλάσσεται στους 2-8°C (35,6-46,4°F)
	- Prodotto da - Fabriqué par - Hergestellt von - Fabricado por	- Manufactured by - Fabricado por - Κατασκευάζεται από
CO-CAL	- Calibratore cut-off - Etalon Seuil - Grenzwert Kalibrator - Calibrador de cut-off	- Cut off Calibrator - Calibrador de cut-off - Οριακός ορός Αντιδραστήριο βαθμονόμησης
CON +	- Controllo positivo - Contrôle Positif - Positiv Kontrolle - Controllo positivo	- Positive Control - Control Positivo - Θετικός ορός ελέγχου
CON -	- Controllo negativo - Contrôle Négatif - Negativ Kontrolle - Controllo negativo	- Negative Control - Control Negativo - Αρνητικός ορός ελέγχου
CAL	- Calibratore - Etalon - Kalibrator - Calibrador	- Calibrator - Calibrador - Αντιδραστήριο βαθμονόμησης
CONJ	- Coniugato - Conjugé - Konjugat - Conjugado	- Conjugate - Conjugado - Σύζευγμα
MP	- Micropietra rivestita - Microplaque sensibilisée - Beschichtete Mikrotiterplatte - Microplaca revestida	- Coated microtiter plate - Microplaca sensibilizada - Επικαλυμμένη μικροπλάκα
WASHB 50x	- Tampone di lavaggio - Tampon de Lavage - Waschpuffer - Solução de lavagem	- Wash buffer - Solución de lavado - Ρυθμιστικό διάλυμα πλύσης
SUB	- Tampone substrato - Substrat - Substratpuffer - Substrato	- Substrate buffer - Tampón sustrato - Ρυθμιστικό διάλυμα υποστρώματος
STOP	- Reagente bloccante - Solution d'Arrêt - Stopreagenz - Solução de paragem	- Stop solution - Solución de parada - Αντιδραστήριο διακοπής αντίδρασης
SB 5x	- Tampone campione - Tampon Echantillons - Probenpuffer - Diluente de amostra	- Sample buffer - Tampón Muestras - Ρυθμιστικό διάλυμα δειγμάτων