
AESKU.DIAGNOSTICS

THE DIAGNOSTIC TOOL THAT WORKS

AESKULISA®

THE DIAGNOSTIC TOOL THAT WORKS

INSTRUCTION MANUAL

AESKULISA® ASCA-A

Ref 3507

Product Ref.	3507
Product Desc.	ASCA-A
Manual Rev. No.	005a: 2025-09-04

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Indice

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12 Bibliografia

1. Seibold F, Stich O, Hufnagel R, Kamil S, Scheurlen M (2001).

Anti-Saccharomyces cerevisiae antibodies in inflammatory bowel disease: a family study.

Scand J Gastroenterol 36: 196-201.

2. Sendid B, Colombel JF, Jacquinot PM et al. (1996).

Specific antibody response to pligomannosidic epitopes in Crohn`s disease.

Clin Diagn Lab Immunol 3: 219-226.

3. Giaffer MH, Clark A, Holdsworth CD (1992).

Antibodies to Saccharomyces cerevisiae in patients with Crohn`s disease and their possible pathogenic importance.

Gut 33: 1071-1075.

4. Bernstein CN, Blanchard JC, Rawsthorne P and Wadja A (1999).

Epidemiology of Crohn`s disease in a central canadian province: a population-based study.

Am. J. Epidemiol. 149(10): 916-924.

5. Glas J, Torok HP, Vilsmaier F, Herbinger KH, Hoelscher M, Folwaczny C (2002)

Anti-saccharomyces cerevisiae antibodies in patients with inflammatory bowel disease and their first-degree relatives: potential clinical value.

Digestion. 66(3):173-177.

- 6. Peters M, Joossens S, Vermeire S, Vlietinck R, Bossuyt X, Rutgeerts P (2001)**
Diagnostic value of Anti-Saccharomyces Cerevisiae and antineutrophil cytoplasmatic autoantibodies in inflammatory bowel disease.
 Am. J. Gastroenterol. 96: 730-734

- 7. Sandborn WJ. (2004)**
Serologic markers in inflammatory bowel disease: state of the art.
 Rev. Gastroenterol. Disord. 4(4):167-74

- 8. Rutgeerts P and Vermeire S (1998)**
Clinical value of the detection of antibodies in the serum for diagnosis and treatment of inflammatory bowel disease.
 Gastroenterol. 115: 1006-1022

- 9. Klebl FH, Bataille F, Hofstädter F, Herfarth H, Schölmerich J, Rogler G (2004)**
Optimising the diagnostic value of Anti-Saccharomyces Cerevisiae-antibodies (ASCA) in Crohn's disease.
 Int. J. Colorectal Dis. 19: 319-24

- 10. Lawrance IC, Murray K, Hall A, Sung JJ, Leong R. (2004)**
A prospective comparative study of ASCA and pANCA in Chinese and Caucasian IBD patients.
 Am. J. Gastroenterol. 99: 2186-94

- 11. Damoiseaux JG, Bouten B, Linders AM, Austen J, Roozendaal C, Russel MG, Forget PP, Tervaert JW. (2002)**
Diagnostic value of anti-Saccharomyces cerevisiae and antineutrophil cytoplasmic antibodies for inflammatory bowel disease: high prevalence in patients with celiac disease.
 J. Clin. Immunol. 22: 281-87

- 12. Muratori P, Muratori L, Guidi M, Maccariello S, Pappas G, Ferrari R, Gionchetti P, Campieri M, Bianchi FB. (2003)**
Anti-Saccharomyces cerevisiae antibodies (ASCA) and autoimmune liver diseases.
 Clin. Exp. Immunol. 132: 473 – 476

- 13. Sutton CL, Yang H, Li Z, Rotter JI, Targan SR, and Braun J. (2003)**
Familial expression of anti-Saccharomyces cerevisiae mannan antibodies in affected and unaffected relatives of patients with Crohn's disease.
 Gut 46: 58 - 63.

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- 15. CLSI Guideline GP44-A4:** Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests

13 Simboli normativi

IVD	- Diagnosi in vitro - Pour diagnostic in vitro - In Vitro Diagnostikum - Para uso Diagnóstico in vitro	- For in vitro diagnostic use - Para uso diagnóstico in vitro - In Vitro Διαγνωστικό μέσο
REF	- Numero d'ordine - Référence Catalogue - Bestellnummer - Número de catálogo	- Catalogue number - Numéro de catálogo - Αριθμός παραγγελίας
LOT	- Descrizione lotto - Lot - Chargen Bezeichnung - Lote	- Lot - Lote - Χαρακτηρισμός παρτίδας
UDI	- Identificatore univoco del dispositivo - Identifiant unique de l'appareil - eindeutige Produktidentifizierung - Identificador único do dispositivo	- Unique Device Identifier - Identificador único del dispositivo - Μοναδικό αναγνωριστικό συσκευής
CE	- Conformità europea - Déclaration CE de Conformité - Europäische Konformität - Declaração CE de Conformidade	- EC Declaration of Conformity - Declaración CE de Conformidad - Ευρωπαϊκή συμφωνία
	- 96 determinazioni - 96 tests - 96 Bestimmungen - 96 Testes	- 96 tests - 96 pruebas - 96 προσδιορισμοί
	- Rispettare le istruzioni elettroniche per l'uso - Voir les instructions d'utilisation électronique - Elektronische Gebrauchsanweisung beachten - Seguir as instruções electrónicas de utilização	- See electronic instructions for use - Siga las instrucciones electrónicas de uso - Ακολουθήστε τις ηλεκτρονικές οδηγίες χρήσης
	- Da utilizzarsi entro - Utilise avant le - Verwendbar bis - Utilizar antes de	- Use by - Utilizar antes de - Χρήση μέχρι
	- Conservare a 2-8°C (35,6-46,4°F) - Conserver à 2-8°C (35,6-46,4°F) - Lagerung bei 2-8°C (35,6-46,4°F) - Conservar entre 2-8°C (35,6-46,4°F)	- Store at 2-8°C (35,6-46,4°F) - Conservar a 2-8°C (35,6-46,4°F) - Φυλάσσεται στους 2-8°C (35,6-46,4°F)
	- Prodotto da - Fabriqué par - Hergestellt von - Fabricado por	- Manufactured by - Fabricado por - Κατασκευάζεται από
CO-CAL	- Calibratore cut-off - Etalon Seuil - Grenzwert Kalibrator - Calibrador de cut-off	- Cut off Calibrator - Calibrador de cut-off - Οριακός ορός Αντιδραστήριο βαθμονόμησης
CON +	- Controllo positivo - Contrôle Positif - Positiv Kontrolle - Controllo positivo	- Positive Control - Control Positivo - Θετικός ορός ελέγχου
CON -	- Controllo negativo - Contrôle Négatif - Negativ Kontrolle - Controllo negativo	- Negative Control - Control Negativo - Αρνητικός ορός ελέγχου
CAL	- Calibratore - Etalon - Kalibrator - Calibrador	- Calibrator - Calibrador - Αντιδραστήριο βαθμονόμησης
CONJ	- Coniugato - Conjugué - Konjugat - Conjugado	- Conjugate - Conjugado - Σύζευγμα
MP	- Micropiastra rivestita - Microplaque sensibilisée - Beschichtete Mikrotiterplatte - Microplaca revestida	- Coated microtiter plate - Microplaca sensibilizada - Επικαλυμμένη μικροπλάκα
WASHB 50x	- Tampone di lavaggio - Tampon de Lavage - Waschpuffer - Solução de lavagem	- Wash buffer - Solución de lavado - Ρυθμιστικό διάλυμα πλύσης
SUB	- Tampone substrato - Substrat - Substratpuffer - Substrato	- Substrate buffer - Tampón sustrato - Ρυθμιστικό διάλυμα υποστρώματος
STOP	- Reagente bloccante - Solution d'Arrêt - Stopreagenz - Solução de paragem	- Stop solution - Solución de parada - Αντιδραστήριο διακοπής αντίδρασης
SB 5x	- Tampone campione - Tampon Echantillons - Probenpuffer - Diluente de amostra	- Sample buffer - Tampón Muestras - Ρυθμιστικό διάλυμα δειγμάτων