




AESKU.DIAGNOSTICS

THE DIAGNOSTIC TOOL THAT WORKS



AESKUBLOTS[®]

THE DIAGNOSTIC TOOL THAT WORKS

INSTRUCTION MANUAL

AESKUBLOTS[®] Gastro Pro

Ref 4005





Product Ref.	4005
Product Desc.	Gastro Pro
Manual Rev. No.	011a: 2025-04-02

Instruction Manual

Table of Contents

1	Intended Use.....	1
2	Clinical Application and Principle of the Test.....	1
3	Kit Contents.....	2
4	Storage and Shelf Life.....	3
5	Precautions of Use and General Introductions.....	3
6	Sample Collection, Handling and Storage.....	4
7	Assay Procedure.....	5
8	Qualitative Interpretation.....	8
9	Technical Data.....	10
10	Performance Data.....	10
11	Disposal.....	16
12	Literature.....	17
13	Regulatory Symbols.....	18



Aesku.Diagnostics GmbH & Co. KG
Mikroforum Ring 2
55234 Wendelsheim
Germany
Phone: +49 6734 9622-0
Fax: +49 6734 9622-2222
Website: www.aesku.com
Mail: info@aesku.com



Product Ref.	4005
Product Desc.	Gastro Pro
Manual Rev. No.	011a: 2025-04-02

1 Intended Use

AESKUBLOTS® Gastro Pro is a membrane-based enzyme immunoassay for the qualitative detection of IgG and IgA antibodies against gliadin, tTg neo-epitope, mannan (ASCA), parietal cell antigen and intrinsic factor in human serum. Antigens are located as parallel lines at exactly defined positions on a nitrocellulose membrane.

The assay is a tool in diagnosis of celiac disease, pernicious anaemia and inflammatory bowel diseases.

This test-kit is for professional use only in laboratories.

2 Clinical Application and Principle of the Test

Each section of the gastrointestinal tract may be affected by autoimmune gastrointestinal diseases. The diseases are often diagnosed years after the first onset of symptoms and in many cases, they have a severe course.

Celiac patients often have an IgA deficiency. In order to avoid false-negative results, this test detects IgA and IgG antibodies.

Antibodies against:

- gliadin are typical for celiac disease. The IgA type is essentially specific for celiac disease. Antibodies of the IgG type occur in 40-50 % of patients with Crohn's disease and also in 10-20 % of patients with ulcerative colitis.
- the neoepitope of tTg (neo-tTg; tissue transglutaminase cross-linked to gliadin specific peptides) constitutes a reliable marker of celiac disease and dermatitis herpetiformis (Dühring's disease). Due to the structural similarity with the physiological epitopes, antibodies against neo-tTg exhibit greater sensitivity (98-100 %) and specificity (93-96 %) than anti-tTg antibodies.
- Mannan (ASCA) have a specificity of 97 % for Crohn's disease. They are essential for the differential diagnosis of Crohn's disease and ulcerative colitis. Up to 75 % of Crohn's disease patients show, in contrast to patients with ulcerative colitis, an increased antibody level.
- parietal cells can be detected by immunofluorescence test in 80-90 % of patients with pernicious anaemia but can also be found in 2-5 % of healthy individuals. They also occur in patients with autoimmune endocrinal diseases and chronic atrophic gastritis type A.
- intrinsic factor show for pernicious anaemia a sensitivity of 50-70 % with a specificity of 100 % in a population of healthy blood donors. Also, these antibodies are detected in patients with autoimmune thyroid diseases and chronic atrophic gastritis type A.



Product Ref.	4005
Product Desc.	Gastro Pro
Manual Rev. No.	011a: 2025-04-02

Principle of the test

The antigens are applied as lines on a nitrocellulose membrane. The membrane is blocked to prevent unspecific reactions. Membrane-strips with specific antigens at exactly defined positions are incubated in serum samples diluted 1:101. Patient's antibodies, if present in the specimen, bind to the antigen. The unbound fraction is washed off in the following step. Afterwards, anti-human immunoglobulins conjugated to horseradish peroxidase (conjugate) are incubated and react with the antigen-antibody complex of the samples. Unbound conjugate is washed off in the following step. After the addition of the TMB-substrate it is converted by an enzymatic reaction to a blue precipitate. The reaction is stopped by distilled water.

3 Kit Contents

TO BE RECONSTITUTED				
Item	Quantity	Cap color	Solution color	Description / Contents
Blocking Reagent	3 x for 10 ml Concentrate each	white	N/A	Non-fat dry milk powder for preparation of 3 x 10 ml sample buffer
Wash Buffer (20x)	1 x 50 ml	white	colorless	20x concentrated for preparation of 1 L Tris buffer, pH 6.9 ± 0.2
READY TO USE				
Item	Quantity	Cap color	Solution color	Description / Contents
Conjugate, IgG	1 x 10 ml	blue	colorless	Anti-human immunoglobulin G (IgG) conjugated to horseradish peroxidase
Conjugate, IgA	1 x 10 ml	red	colorless	Anti-human immunoglobulin A (IgA) conjugated to horseradish peroxidase
TMB Substrate	1 x 10 ml	black	colorless	Stabilized TMB/H ₂ O ₂
Membrane strips	24 strips	color coding: black	N/A	Coated antigens see Intended use
tweezers, reference template, scoring sheet, adhesive strip (double-sides, white)	1 pcs. each	N/A	N/A	N/A
incubation tray	3 pcs.	N/A	N/A	N/A
Labels for sample buffer	3 pcs.	N/A	N/A	N/A
MATERIALS REQUIRED, BUT NOT PROVIDED				
rocking platform, cylinder 1000 ml, pipette or cylinder for 10 ml, precision pipettes (10, 1000 µl), absorbent or filter paper. Our tests are designed to be used with purified water according to the definition of the United States Pharmacopeia (USP 26 - NF 21) and the European Pharmacopeia (Eur.Ph. 4th ed.).				



Product Ref.	4005
Product Desc.	Gastro Pro
Manual Rev. No.	011a: 2025-04-02

4 Storage and Shelf Life

Store all reagents and membrane-strips at 2-8 °C/35.6-46.4 °F in their original containers. Once prepared, reconstituted wash buffer as well as opened strips, conjugate and TMB are stable at 2-8 °C/35.6-46.4 °F for six weeks. Reconstituted blocking reagent is stable at 2-8 °C/35.6-46.4 °F for 3 weeks. Reagents and strips shall be used within the expiry date indicated on each respective component. Don't use components after the expiry dates. Avoid intense exposure of TMB solution to the light.

5 Precautions of Use and General Introductions

5.1 Health hazard data

This product is for IN VITRO DIAGNOSTIC use only. Thus, only staff trained and specially advised in methods of in vitro diagnostics may perform the kit.

All kit components are classified according to Regulation (EC) No. 1272/2008 [CLP]. See the Materials and Safety Document (MSDS) for more information on ingredients.

Substances listed on the so-called "Candidate List of Substances of very High Concern (SVHCV) for authorization" of the European Chemicals Agency (ECHA) are not intentional components of this product. It is therefore not to be expected that these substances are contained in amounts $\geq 0.1\%$ in the product.

Reagents should be stored safely and be inaccessible to children.


In particular, the mixture does not contain any substances in concentrations $\geq 0.1\%$ to be classified as PBT or vPvB.

Patient samples should be considered potentially infectious and handled according to national laws. Patient samples and other potentially infectious material should be decontaminated after the test run.

5.2 General directions for use

To differentiate between the various **AESKUBLOTS**[®]-tests available, a color coding is applied above the reference line of the strips:

Color coding	AESKUBLOTS [®]
red	ANA-17 comp
orange	ANA-17 Pro
brown	Liver Pro
purple	Vasculitis Pro
black	Gastro Pro
green	Borrelia-G and Borrelia-M

 AESKU DIAGNOSTICS THE DIAGNOSTIC TOOL THAT WORKS	Product Ref.	4005
	Product Desc.	Gastro Pro
	Manual Rev. No.	011a: 2025-04-02

In case the product information, including the labeling, is incorrect please contact the manufacturer or the supplier of the test kit.

Blocking Reagent and wash buffer may be interchanged between lots and test kits. All other components are specific for each test kit and are not to be interchanged. Don't exchange substances between autoimmunity and borrelia diagnostic tests!

For handling of conjugate do not use polystyrene vessels.

Allow all components to reach room temperature (20-32 °C/68-89.6 °F) before use, mix well and follow the recommended incubation scheme for an optimum performance of the test.

Never expose components to higher temperature than 37 °C/ 98.6 °F.

Always pipette substrate solution with brand new tips only. Protect this reagent from light. Never pipette conjugate with tips priorly used with other reagents.

The intensity of the band color does not necessarily correlate with antibody titers obtained by other reference methodologies.

Samples from apparent normal blood donors may contain autoantibodies.

If the patient sample contains elevated levels of immune complexes or other immunoglobulin aggregates, false positive results by non-specific binding cannot be ruled out.

A definite clinical diagnosis should not be based on the results of the performed test only but should be made by the physician after all clinical and laboratory findings have been evaluated. The diagnosis is to be verified using different diagnostic methods.

6 Sample Collection, Handling and Storage

Use preferentially freshly collected serum samples. Blood withdrawal must follow national requirements. Do not use icteric, lipemic, hemolyzed or bacterially contaminated samples. Sera with particles should be cleared by low speed centrifugation (< 1000 x g). Blood samples should be collected in clean, dry and empty tubes.

After separation, the serum samples should be used during the first 8 h. Alternatively, the samples should be stored in tightly closed vials at 2-8 °C/35.6-46.4 °F for up to 48 h, or frozen at -20 °C/-4 °F for longer periods. Avoid freezing and thawing samples repeatedly. (Thomas: Labor und Diagnose; CLSI Guideline GP44-A4) Do not use thermally inactivated samples (56 °C/132.8 °F).



Product Ref.	4005
Product Desc.	Gastro Pro
Manual Rev. No.	011a: 2025-04-02

7 Assay Procedure

7.1 Preparations prior to starting

Confirm that no salt crystals have been formed in the concentrate. If this happened, dissolve the crystals by slightly warming, room temperature should be enough, the concentrate.

Dilute concentrated wash buffer 1:20 with distilled water (e.g. 50 ml plus 950 ml).

For preparation of sample buffer: add 10 ml wash buffer to one bottle Blocking Reagent and mix well.

7.2 Test Steps

Important notes:

Follow exactly this protocol. Make sure that the two components mentioned in the protocol are added to the tray in step 2, 6, 9.








Do not let strip dry out during incubation steps.

Do not touch strip with fingers, use tweezers.




Remove diluted samples completely after incubation of strip to avoid carry over.

Continuously shake strip during incubation steps.

Give sample buffer, conjugate and substrate together with the wash buffer to one side of the incubation tray. Do not allow to flow over the strip.

Step	Description
1.	Ensure the preparations, from step 7.1 above, have been carried out prior to test begin.
2.	 <p>Put strip in correct orientation into incubation tray (reference line and color coding upwards). Put 700 μl wash buffer and 300 μl sample buffer in the incubation tray. Moisten strip with the solution and incubate for 5 minutes with agitation.</p>
CONTROLS & SAMPLES	
3.	 <p>Samples</p> <p>Pipette 10 μl serum sample into the designated incubation trays with sample buffer.</p>
4.	 <p>Incubate for 30 minutes at 20-32 °C/68-89.6 °F with agitation. After that remove sample completely.</p>
5.	 <p>Wash 3 times for 5 minutes with 1.5 ml wash buffer by agitation. Remove wash buffer after every washing step.</p> <p style="text-align: center;">Σ 3x</p>
CONJUGATE	
6.	 <p>Pipette 700 μl wash buffer and 300 μl conjugate into each incubation tray with strip.</p>
7.	 <p>Incubate for 30 minutes at 20-32 °C/68-89.6 °F with agitation. Remove conjugate.</p>
8.	 <p>Wash 3 times for 5 minutes with 1.5 ml wash buffer by agitation. Remove wash buffer after every washing step.</p> <p style="text-align: center;">Σ 3x</p>

SUBSTRATE

9.		Pipette 700 µl dH ₂ O and 300 µl substrate into each incubation tray with strip.
10.		Incubate for 15 minutes at 20-32 °C/68-89.6 °F with agitation, protected from intense light. Remove substrate.
STOP		
11.		Pipette 2 ml dH ₂ O into each incubation tray with strip. Incubate 1 minute with agitation. Remove dH ₂ O. Repeat this step one time.
12.	Remove strip of the incubation tray. Dry strip between filter paper	
13.	Analyze results within 24 h.	

AESKUBLOTS[®] Gastro Pro is also intended to be automatically processed and evaluated on the **HELIA[®] Automated Blot System**.

Reagent preparation for **HELIA[®] Automated Blot System**: Dilute one part wash buffer concentrate (WASH) with 19 parts ultrapure water (e.g. 50 ml wash buffer concentrate and 950 ml ultrapure water) to obtain a ready-to-use wash buffer. All other reagents are ready to use when processed in **HELIA[®] Automated Blot System**. For detailed handling of the test on **HELIA[®] Automated Blot System** refer to the instruction manual of the **HELIA[®] Automated Blot System**.



Product Ref.	4005
Product Desc.	Gastro Pro
Manual Rev. No.	011a: 2025-04-02

8 Qualitative Interpretation

In case that the values of the controls do not meet the criteria, the test is invalid and has to be repeated. We recommend retesting samples that are equivocal.

The following technical issues should as well be checked: expiry date of (prepared) reagents, storage conditions, pipettes, equipment, incubation conditions and washing methods.

If the tested samples show aberrant values or any kind of deviation or if the validation criteria are not met because of reasons outside the operator's responsibility, please contact the manufacturer or the supplier of the test kit.

Medical laboratories might perform an in-house quality control by using their own controls and/or internal pooled sera, as stated in national regulations.

8.1 Manual Analysis

Test results can be considered valid, if:


- Functional control is visible
- Cut-off control is visible
- Color intensity of Cut-off control is weaker than color intensity of functional control

Fix dried strip onto scoring sheet aligned with reference line. Align reference template with the strip reference line. Interpret results only in reference to Cut-off control of each strip.

Each test kit contains a color copy with all bands provable in the test.

The analysis is carried out by means of comparing the color intensities of the bands with color intensity of the Cut-off control. The test is equivocal if the intensities do not significant differ. If the color is more intense, the test result is positive, if the color intensity is weaker, the test is negative.

The results can be recorded on the scoring sheet.

 AESKU.DIAGNOSTICS <small>THE DIAGNOSTIC TOOL THAT WORKS</small>	Product Ref.	4005
	Product Desc.	Gastro Pro
	Manual Rev. No.	011a: 2025-04-02

8.2 Software-supported evaluation

The analysis of the strips can be carried out by means of using AESKU.SCAN Software. Please refer to the instructions for use of AESKU.SCAN.

Test results can be considered valid, if:

- Functional control is visible
- Cut-off control is visible
- Color intensity of Cut-off control is weaker than color intensity of functional control

AESKU.SCAN 2.0:

Fix dried strip onto scoring sheet (printable) aligned with reference line. Align reference template with the strip reference line.

Evaluate strips according to the instructions for use of AESKU.SCAN 2.0 software.

Qualitative result analysis is carried out by means of comparing the color intensities of the individual antigens with the color intensity of the Cut-off control.

AESKU.SCAN 3.0:

Put strips within the incubation tray into the reader.

Evaluate strips according to the instructions for use of AESKU.SCAN 3.0 software.

Qualitative result analysis is carried out by means of comparing the color intensities of the individual antigens with the color intensity of the Cut-off control.

HELIA® Automated Blot System:

Using a **HELIA® Automated Blot System**, the results are analyzed automatically. The results can be determined in Index-values.

The following interpretation according to the signal intensity is suggested:

Result Interpretation	Symbol	Index	Color
Negative	-	0.0 - <0.8	Colorless
Equivocal	+/-	≥0.8 - <1.15	Blue
Weak positive	+	≥1.15 - <2.5	Yellow
Positive	++	≥2.5 - <4.0	Red
Strong positive	+++	≥ 4.0	Dark red



Product Ref.	4005
Product Desc.	Gastro Pro
Manual Rev. No.	011a: 2025-04-02

9 Technical Data

Sample material:	serum
Sample volume:	10 µl of sample
Total incubation time:	112 minutes at 20-32 °C/68-89.6 °F
Storage:	at 2-8 °C/35.6-46.4 °F; use original vials only.
Number of determinations:	24 tests

10 Performance Data

10.1 Normal Range Study

Expected values for **AESKUBLOTS® Gastro Pro** were analyzed with a panel of 120 sera from healthy donors.


All assays were performed fully automated according to the current IFUs.

AESKUBLOTS® Gastro Pro IgA

AESKUBLOTS® Gastro Pro - REF 4005- Normal Range- Report										
Antigen	Number of samples	Positive samples		Equivocal samples		Negative samples		Min [Index]	Max [Index]	Mean [Index]
		n	[%]	n	[%]	n	[%]			
Gliadin	120	0	0.00	9	7.50	111	92.50	0.00	1.08	0.27
Neo-tTG	120	1	0.83	0	0.00	119	99.17	0.00	1.30	0.03
ASCA	120	0	0.00	11	9.17	109	90.83	0.00	1.10	0.28
PCA	120	0	0.00	10	8.33	110	91.67	0.00	1.08	0.24
IF	120	0	0.00	0	0.00	120	100.00	0.00	0.70	0.04

AESKUBLOTS® Gastro Pro IgG

AESKUBLOTS® Gastro Pro - REF 4005- Normal Range- Report										
Antigen	Number of samples	Positive samples		Equivocal samples		Negative samples		Min [Index]	Max [Index]	Mean [Index]
		n	[%]	n	[%]	n	[%]			
Gliadin	120	0	0.00	8	6.67	112	93.33	0.00	1.10	0.21
Neo-tTG	120	1	0.83	0	0.00	120	100.00	0.00	0.21	0.00
ASCA	120	0	0.00	5	4.17	115	95.83	0.00	0.98	0.30
PCA	120	0	0.00	4	3.33	116	96.67	0.00	1.03	0.14
IF	120	0	0.00	0	0.00	120	100.00	0.00	0.27	0.01

 AESKU.DIAGNOSTICS THE DIAGNOSTIC TOOL THAT WORKS	Product Ref.	4005
	Product Desc.	Gastro Pro
	Manual Rev. No.	011a: 2025-04-02

With **AESKUBLOTS® Gastro Pro IgA** there was one positive result for anti-Neo-tTG. The most equivocal results were detected for anti-Gliadin in both variants. The lowest Index-value measured was 0.0, the highest Index-value measured was 1.3 (Neo-tTG) for **AESKUBLOTS® Gastro Pro IgA**.

The low number of positive samples we found with **AESKUBLOTS® Gastro Pro** in the tested healthy population correlates well with the numbers reported in the literature.

We also recommend in the IFU in chapter 8: “Qualitative interpretation” that each laboratory should establish its own normal range.

10.2 Precision

Intra-Assay Precision

To determine the within-lab / intra-assay variances four different sera (S1-S4) have been tested for 20 times on **AESKUBLOTS® Gastro Pro**. The positive- and negative agreement were calculated from these test values (n=20).

AESKUBLOTS® Gastro Pro IgA

Calculation of agreements	
Percent positive agreement [%]	100.0
Percent negative agreement [%]	100.0
Percent overall agreement [%]	100.0

AESKUBLOTS® Gastro Pro IgG

Calculation of agreements	
Percent positive agreement [%]	100.0
Percent negative agreement [%]	100.0
Percent overall agreement [%]	100.0

Inter-Assay Precision

To determine the inter-assay / day to day variances four different sera (S1-S4) over the whole range have been tested 8 times per run, in a total of 5 runs on 5 different days (n=40). For day 1 the results from the within-lab / intra-assay variances are used (n=20), the positive- and negative agreement were calculated from these test values.

AESKUBLOTS® Gastro Pro IgA

Calculation of agreements	
Percent positive agreement [%]	100.0
Percent negative agreement [%]	100.0
Percent overall agreement [%]	100.0



Product Ref.	4005
Product Desc.	Gastro Pro
Manual Rev. No.	011a: 2025-04-02

AESKUBLOTS® Gastro Pro IgG

Calculation of agreements	
Percent positive agreement [%]	100.0
Percent negative agreement [%]	100.0
Percent overall agreement [%]	100.0

Lot to Lot Precision

To determine the Lot-to-Lot Precision, four different sera (S1-S4) over the whole range have been tested 8 times on 3 Lots (n=24). Lot 1 is taken from Run 1 of the Inter Assay data. The positive and negative agreement were calculated across the test values of the three different Lots.

AESKUBLOTS® Gastro Pro IgA

Calculation of agreements	
Percent positive agreement [%]	100.0
Percent negative agreement [%]	100.0
Percent overall agreement [%]	100.0

AESKUBLOTS® Gastro Pro IgG

Calculation of agreements	
Percent positive agreement [%]	100.0
Percent negative agreement [%]	100.0
Percent overall agreement [%]	100.0

10.3 Diagnostic Sensitivity and Specificity

A clinical study with different serum samples from patients with celiac disease, Dermatitis herpetiformis, Morbus Crohn, Ulcerative Colitis, SLE, RA, Diabetes type 1 and healthy controls were analyzed with the **AESKUBLOTS® Gastro Pro**. Furthermore, these samples were measured with a competitor western blot. Also, a pediatric sample cohort with clear diagnosis and medical history of 87 samples was included.

All assays were performed according to the current IFUs. After that the results were compared to determine the diagnostic sensitivity and specificity of the assay.

Anti-Gliadin

Gliadin are typical for celiac disease. The IgA type is essentially specific for celiac disease. Antibodies of the IgG type occur in 40-50 % of patients with Crohn's disease and also in 10-20 % of patients with ulcerative colitis.



Product Ref.	4005
Product Desc.	Gastro Pro
Manual Rev. No.	011a: 2025-04-02

AESKUBLOTS® Gastro Pro IgA and Aeskublots® Gastro Pro IgG

Celiac samples were analyzed, however, no sensitivity for gliadin can be calculated as the reference kit does not have gliadin for IgA and IgG variant. The prevalence of celiac disease-specific autoantibodies has been studied in patients with SLE, RA, and diabetes type 1.

	n (65)	AESKUBLOTS® Gastro Pro IgA anti-gliadin positive	AESKUBLOTS® Gastro Pro IgG anti-gliadin positive
SLE	28	7.1 %	10.7 %
RA	11	0.0 %	0.0 %
Diabetes type 1	26	0.0 %	34.6 %

Gliadin antibodies: analysis of healthy controls	AESKUBLOTS® Gastro Pro IgA		
	Positive	Equivocal	Negative
	0	9	111
Gliadin antibodies: analysis of healthy controls	AESKUBLOTS® Gastro Pro IgG		
	Positive	Equivocal	Negative
	0	8	112

In the analysis of 120 antibody-negative healthy patient samples, no positive samples were measured, and only 9 and 8 samples, respectively, showed equivocal results. It is known from the literature that gliadin antibodies can also occur in healthy patients or in patients with other diseases.

Anti-tTG neo

AESKUBLOTS® Gastro Pro IgA

Neo tTG 2x2	Diagnose	
	Celiac Diseases	Controls
Positive	28	1
Negative	18	40

Diagnostic performance tTG-neo IgA	CI
Sensitivity	60.9 45.4 % to 74.9 %
Specificity	97.6 87.1 % to 99.9 %
AUC	0.792 0.692 to 0.872

In a pediatric sample cohort of 87 samples including Celiac Disease, Morbus Chron, Ulcerative Colitis and healthy blood donors a diagnostic sensitivity for tTG-neo IgA of 60.9 % and a diagnostic specificity of 97.6 % could be calculated.



Product Ref.	4005
Product Desc.	Gastro Pro
Manual Rev. No.	011a: 2025-04-02

AESKUBLOTS® Gastro Pro IgG

Neo tTG 2x2	Diagnose	
	Celiac Diseases	Controls
Positive	38	3
Negative	10	33

Diagnostic performance tTG-neo IgA		CI
Sensitivity	79.2	65.1 % to 89.5 %
Specificity	91.7	77.5 % to 98.2 %
AUC	0.854	0.760 to 0.922

In a pediatric sample cohort of 84 samples including Celiac Disease, Morbus Chron, Ulcerative Colitis and healthy blood donors a diagnostic sensitivity for tTG-neo IgG of 79.2 % and a diagnostic specificity of 91.7 % could be calculated.

Anti-ASCA

AESKUBLOTS® Gastro Pro IgA

The analysis of 24 sera from patients with Crohn`s disease (anti-ASCA positive n = 15) showed a sensitivity of 100 % for the detection of autoantibodies against ASCA with respect to the reference method Western blot of a competitor. In comparison to the competitor assay a specificity of 100 % could be measured.

ASCA antibodies: analysis of anti-ASCA in Crohn`s disease patients		Competitor assay		
		Positive	Equivocal	Negative
AESKUBLOTS® Gastro Pro IgA	Positive	15	1	0
	Equivocal	2	1	0
	Negative	0	1	4

The prevalence of ASCA autoantibodies has been studied in patients with Crohn`s disease, SLE, RA, diabetes type 1 and healthy controls.

	n (124)	AESKUBLOTS® Gastro Pro IgA anti-ASCA positive
Crohn's disease	24	70.0 %
SLE	28	0.0 %
RA	11	0.0 %
Diabetes type 1	26	25.0 %
Healthy controls	35	0.0 %



Product Ref.	4005
Product Desc.	Gastro Pro
Manual Rev. No.	011a: 2025-04-02

AESKUBLOTS® Gastro Pro IgG

The analysis of 20 sera from patients with Crohn`s disease (anti-ASCA positive n = 9) showed a sensitivity of 100 % for the detection of autoantibodies against ASCA with respect to the reference method Western blot of a competitor. In comparison to the competitor assay a specificity of 100 % could be measured.

ASCA antibodies: analysis of anti-ASCA in Crohn`s disease patients		Competitor assay		
		Positive	Positive	Positive
AESKUBLOTS® Gastro Pro IgG	Positive	9	1	0
	Equivocal	0	3	3
	Negative	0	1	3

The prevalence of ASCA autoantibodies has been studied in patients with Crohn`s disease, SLE, RA, diabetes type 1 and healthy controls.

	n (136)	AESKUBLOTS® Gastro Pro IgG anti-ASCA positive
Crohn`s disease	20	50.0 %
SLE	28	7.1 %
RA	11	9.1 %
Diabetes type 1	18	5.6 %
Healthy controls	59	0.0 %

Anti-PCA

Parietal cells can be detected by immunofluorescence test in 80-90 % of patients with pernicious anaemia but can also be found in 2-5 % of healthy individuals. They also occur in patients with autoimmune endocrinal diseases and chronic atrophic gastritis type A.

In the sample panel used, there were no samples of pernicious anaemia, in which antibodies against PCA typically occur. Therefore, no sensitivity can be given here for **AESKUBLOTS® Gastro Pro** IgA and IgG. In sera from patients with SLE, RA, diabetes type 1 and healthy controls the prevalence for anti-PCA was studied.

	n (185)	AESKUBLOTS® Gastro Pro IgA anti-PCA positive	AESKUBLOTS® Gastro Pro IgG anti-PCA positive
SLE	28	10.7 %	10.7 %
RA	11	9.1 %	18.2 %
Diabetes type 1	26	26.9 %	23.0 %
Healthy controls	120	0.0 %	0.0 %



Product Ref.	4005
Product Desc.	Gastro Pro
Manual Rev. No.	011a: 2025-04-02

Anti-Intrinsic Factor

Intrinsic factor show for pernicious anaemia a sensitivity of 50-70 % with a specificity of 100 % in a population of healthy blood donors. Also, these antibodies are detected in patients with autoimmune thyroid diseases and chronic atrophic gastritis type A.

In the sample panel used, there were no samples of pernicious anaemia, in which antibodies against Intrinsic Factor typically occur. Therefore, no sensitivity can be given here for **AESKUBLOTS® Gastro Pro** IgA and IgG. In sera from patients with SLE, RA, diabetes type 1 and healthy controls the prevalence for anti-PCA was studied.

	n (185)	AESKUBLOTS® Gastro Pro IgA anti-IF positive	AESKUBLOTS® Gastro Pro IgG anti-IF positive
SLE	28	0.0 %	3.5 %
RA	11	0.0 %	0.0 %
Diabetes type 1	26	0.0 %	0.0 %
healthy controls	120	0.0 %	0.0 %

The sensitivities and specificities as well as the competitor agreements show very good results. In some cases, the sensitivity could not be calculated for individual parameters due to the lack of positives samples in the cohorts. Here further samples are to be investigated in the future.

10.4 Relative Sensitivity and Specificity

In order to determine the positive agreement (relative sensitivity), 30 sera from IIF or ELISA antibody-positive patients were tested in **AESKUBLOTS® Gastro Pro**. For determination of the negative agreement (relative specificity), 100 sera from blood donors were analyzed.

Antigen	positive agreement (relative sensitivity)	negative agreement (relative specificity)
Gliadin	100 %	100 %
Neo-tTg	100 %	100 %
Mannan (ASCA)	100 %	100 %
Parietal cell antigen	100 %	100 %
Intrinsic factor	100 %	100 %

11 Disposal

Please observe the relevant statutory requirements!



Product Ref.	4005
Product Desc.	Gastro Pro
Manual Rev. No.	011a: 2025-04-02

12 Literature

Carmel R (1992). Reassessment of the relative prevalences of antibodies to gastric parietal cell and to intrinsic factor in patients with pernicious anaemia: influence of patient age and race. *Clin Exp Immunol.* 89(1):74-77.

Lerner A (2014). Serological Diagnosis of Celiac Disease –Moving Beyond the Tip of the Iceberg. *International Journal of Celiac Disease.* 2(2): 64-66.

Catassi C, Rätsch I-M, Fabiani E, Rossini M, Coppa GV, Giorgi PL, Bordicchia F, Candela F (1994). Celiac disease in the year 2000: exploring the iceberg. *The Lancet.* 343:200–203.

Holtmeier W, Caspary WF (1998). Antikörperdiagnostik bei Sprue/Zöliakie. *Z. Gastroenterol.* 36: 587–597.

Mäki M, Collin P (1997). Coeliac disease. *Lancet.* 349: 1755–1759.

Oh R, Brown DL (2003). Vitamin B12 deficiency. *Am Fam Physician.* 67: 979–986.

Shan L, Molberg O, Parrot I, Hausch F, Filiz F, Gray GM, Sollid LM, Khosla C (2002). Structural basis for gluten intolerance in Celiac Sprue. *Science* 297:2275–2279.

Toh Ban-Hock, Alderuccio, F. (2004). Pernicious anaemia. *Autoimmunity.* 37: 357–361.

Toh Ban-Hock, van Driel Ian R. (1997). Pernicious anaemia. *NEJM.* 337:1441–1448.



Product Ref.	4005
Product Desc.	Gastro Pro
Manual Rev. No.	011a: 2025-04-02

13 Regulatory Symbols

	" Diagnosi in vitro	" For in vitro diagnostic use
	" Pour diagnostic in vitro	" Para uso diagnóstico in vitro
	" In Vitro Diagnostikum	" In Vitro Διαγνωστικό μέσο
	" Para uso Diagnóstico in vitro	
	" Numero d'ordine	" Catalogue number
	" Référence Catalogue	" Numéro de catálogo
	" Bestellnummer	" Αριθμός παραγγελίας
	" Número de catálogo	
	" Descrizione lotto	" Lot
	" Lot	" Lote
	" Chargen Bezeichnung	" Χαρακτηρισμός παρτίδας
	" Lote	
	" Identificatore univoco del dispositivo	" Unique Device Identifier
	" Identifiant unique de l'appareil	" Identificador único del dispositivo
	" eindeutige Produktidentifizierung	" Μοναδικό αναγνωριστικό συσκευής
	" Identificador único do dispositivo	
	" Conformità europea	" EC Declaration of Conformity
	" Déclaration CE de Conformité	" Declaración CE de Conformidad
	" Europäische Konformität	" Ευρωπαϊκή συμφωνία
	" Declaração CE de Conformidade	
	" 24 determinazioni	" 24 tests
	" 24 tests	" 24 pruebas
	" 24 Bestimmungen	" 24 προσδιορισμοί
	" 24 Testes	
	" Rispettare le istruzioni elettroniche per l'uso	" See electrical instructions for use
	" Voir les instructions d'utilisation électronique	" Siga las instrucciones electrónicas de uso
	" Elektronische Gebrauchsanweisung beachten	" Ακολουθήστε τις ηλεκτρονικές οδηγίες χρήσης
	" Seguir as instruções electrónicas de utilização	
	" Da utilizzarsi entro	" Use by
	" Utilise avant le	" Utilizar antes de
	" Verwendbar bis	" Χρήση μέχρι
	" Utilizar antes de	
	" Conservare a 2-8°C (35.6-46.4°F)	" Store at 2-8°C (35.6-46.4°F)
	" Conserver à 2-8°C (35.6-46.4°F)	" Conservar a 2-8°C (35.6-46.4°F)
	" Lagerung bei 2-8°C (35.6-46.4°F)	" Φυλάσσεται στους 2-8°C (35.6-46.4°F)
	" Conservar entre 2-8°C (35.6-46.4°F)	
	" Prodotto da	" Manufactured by
	" Fabriqué par	" Fabricado por
	" Hergestellt von	" Κατασκευάζεται από
	" Fabricado por	
	" Strip di nitrocellulosa rivestita	" Coated nitrocellulose strip
	" Strip de nitrocellulose couché	" Tira de nitrocelulosa recubierta
	" Nitrozellulosemembran-Streifen mit aufgebracht Antigenen	" Επίστρωση λωρίδα νιτροκυτταρίνης
	" Tira de nitrocelulose revestido	
	" Tampone di lavaggio	" Wash buffer
	" Tampon de Lavage	" Solución de lavado
	" Waschpuffer	" Ρυθμιστικό διάλυμα πλύσης
	" Solução de lavagem	
	" Reagente bloccante	" Blocking Reagent
	" réactif de blocage	" Reactivo bloqueante
	" Blockier-Reagenz	" Αντιδραστήριο αποκλεισμού
	" Bloqueio de reagente	
	" Ricostituire con 10 mL	" Reconstitute with 10 mL
	" reconstituer avec 10 mL	" reconstituir con 10 mL
	" rekonstituieren mit 10 mL	" Ανασύσταση με 10 mL
	" reconstituir com 10 mL	
	" Tampone campione	" Sample buffer
	" Tampon Echantillons	" Tampón Muestras
	" Probenpuffer	" Ρυθμιστικό διάλυμα δειγμάτων
	" Diluente de amostra	
	" Coniugato	" Conjugate
	" Conjugé	" Conjugado
	" Konjugat	" Σύζευγμα
	" Conjugado	
	" Tampone substrato	" Substrate buffer
	" Substrat	" Tampón sustrato
	" Substratpuffer	" Ρυθμιστικό διάλυμα υποστρώματος
	" Substrato	