INTENDED USE

AESKUCARE® Gluten Related Disorders IgA is a near-patient (point of care) in-vitro-immunoassay intended for the qualitative and quantitative simultaneous measurement of human specific IgA antibodies against human tissue transglutaminase (tTG), neo-epitopes of human tissue transglutaminase (tTG neo), microbial transglutaminase (mTG), neo-epitopes of microbial transglutaminase (mTG neo), deamidated gliadin-specific peptides (DGP), gliadin, peptic-tryptic digests of gliadin (Frazer's fraction), human epidermal transglutaminase (TG3) and total IgA in heparinized or Na-EDTA venous or capillary whole blood, plasma or serum. The test results, in conjunction with the patient's medical history and other findings, can be used by professionals to make a nutritional recommendation or support a diagnosis of gluten related disorders (GRD). The test should be performed by qualified personnel who are familiar with in vitro diagnostic methods.

TEST PRINCIPLE

AESKUCARE® Gluten Related Disorders IgA is a modified ELISA (Enzyme Linked Immunosorbent Assay) system. The device contains 4 rows (test areas), where standards, controls and antigens for triplicate measurement are bound to the membrane. Specific IgA (sIgA) antibodies from the sample bind to the corresponding antigens. The bound human IgA is detected by specific anti-human-IgA antibodies. A secondary antibody coupled to alkaline phosphatase is added for signal amplification. After adding the color reagent bound IgA is identified by an indigo colored precipitate. sIgA concentrations can be determined quantitatively by using a separately available reader or qualitatively by the naked eye. The result is read out by comparing the color intensity of the antigen signal (average of triplicate measurement) against the standard signals of the first and fourth test rows (quantitative evaluation) or against the signal intensity of the cut-off calibrator in the first test row (qualitative evaluation). sIgA concentrations are expressed in arbitrary units U/ml (quantitative evaluation) or classified into three levels (qualitative evaluation).

PRECAUTIONS

1. For human in vitro diagnostic use only. 2. Please read the entire contents of these Instructions for Use before performing the test. 3. Do not use reagents beyond their expiration dates. 4. In case of damage to the packaging please inspect the protecting bag of the device carefully for damage. Make sure the reagents are not damaged or open. In case of doubt do not use the test kit to avoid incorrect results or wrong diagnosis. 5. It is not recommended to pool any reagents. 6. The color reagent contains 5-bromo-4-chloro-3'-indolylphosphate and nitro-blue tetrazolium chloride. Antibody solution 1A and sample blocking dilutor contain ProClin below concentration limit (skin sensitization). Avoid contact with skin. Wear suitable gloves. For more information see Material Safety Data Sheet (available on request). 7. All components are single use only. 8. Test results will not be valid if these precautions are not followed. 9. Fluid paths can be blocked by coagulated blood. 10. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

STORAGE

Store in a dark and cool place at 2-8°C/36-46°F. All reagents are ready-to-use, single dose and packaged in screw tubes.

EXPIRATION

The expiration date is printed on labels placed on the test packaging box. The expiration date of the kit is valid for all kit components, even if expiry of single components is different! After expiry all test components have to be discarded.

REAGENTS AND MATERIALS

Materials provided: Instructions for use, evaluation sheet, control certificate (CoA), 2D-barcode control chart for reader, 1x AESKUCARE® Gluten Related Disorders IgA device, 1x reagent rack containing 9 color coded vials, 1x blood sampling pouch containing: 1x disinfection swab **CE** 1x lancet, sterile **CE** 0124, 1x blood sampling capillary 20 µl with lithium heparin **CE**, 1x swab **CE** 1x band-aid **CE**, 1x syringe.

1. DILB – Sample blocking dilutor, 380µl, 1x (lilac lid) 2. WASH - Washing solution, 1ml each, 3x (blue lid) 3. Ab1A – Antibody solution 1A, 800µl, 1x (white lid, translucent tube) (contains <1% polyclonal anti-human-IgA antibody coupled to digoxigenin) 4. Ab2 - Antibody solution 2, 800µl, 1x (yellow lid) (contains <1% mouse monoclonal anti-digoxigenin antibody coupled to alkaline phosphatase) 5. COLR - Color reagent, 800µl each (contains <1% 5-bromo-4-chloro-3'-indolylphosphate/nitro-blue tetrazolium chloride), 2x (white lid, brown tube) 6. STOP - Stop buffer, 1ml, 1x (green lid). **Note:** All reagents are ready for use; store at 2-8°C until expiry date.

Human material used in the AESKUCARE® Gluten Related Disorders IgA device (dried material on the test strips, no infection risk) was tested and found to be negative for HIV, Hepatitis B and Hepatitis C.

Additional materials required but not included in the kit: disposable gloves, timer. For venous whole blood and/or plasma samples: normal venipuncture supplies, heparinized blood collection tubes, laboratory pipette.

SAMPLE MATERIAL

Capillary blood from a finger prick, from earlobes or heels (in case of infants) can be used in this test. The blood sampling capillary provided in the kit should be used to collect the correct volume (20 µl). Heparinized venous blood (20 µl), plasma or serum (8 µl) may also be used as specimen for this test. In case of plasma or serum, laboratory pipettes should be used to dispense the required volume. Please follow the instructions in section "Test Procedure" to collect the sample.

TEST PROCEDURE

Preparation: All reagents and the assay device must be adjusted to room temperature (18-30°C) for at least 30 minutes prior to use. 1. Open the kit and take out the device, the syringe, the reagent rack and the blood sampling pouch. Place them on a flat surface. 2. Open and assemble the reagent rack. 3. Label the device with the patient ID in the grey label area.

Collection of capillary blood sample: Open the blood sampling pouch and take out its contents: alcohol wipe, safety lancet, blood sampling capillary 20µl, gauze and band-aid. Place all items on a clean, flat surface. 1. Massage the fingertip or put it into warm water for a few moments before puncturing so as to optimize blood flow. 2. Clean the site to be pricked with the alcohol wipe. 3. Unlock the safety lancet as shown (fig. step A) and pull off the protective cap. 4. Place the safety lancet on the finger side to be pricked (fig. step B) and press the button firmly. **Note!** Make sure the hand is below the heart level. If the blood flow is initially slow, gently massage the finger and check if the hand is below (!) the level of heart.

step A
step B

5. Collect the drops of blood with the tip of the blood sampling capillary. Hold the blood sampling capillary horizontal or at a slight downward angle and let the blood flow until it is properly filled. Make sure not to cover the (air) hole in the plunger. The blood stop sets the sample volume (20 µl). **CAUTION!** Make sure that you use 20 µl blood for the test, otherwise you will not get valid results. 6. Cover the wound with the band-aid. 7. Immediately open the sample blocking dilutor vial (lilac lid) and put the tip of the blood sampling capillary in a position just above the vial opening. Cover the hole at the end of the plunger with your finger. Gently press down the plunger to add the blood into the sample blocking dilutor vial. Close the sample blocking dilutor vial and mix gently by hand or vortex briefly.

CAUTION! Do not store the blood after adding it to the sample blocking dilutor. Once added to the sample blocking dilutor, all samples have to be tested immediately. If untreated blood is used (i.e. without adding whole blood to sample blocking dilutor) then the test results will not be valid. Do not store the blood within the capillaries.

Collection of venous blood sample, plasma sample, or serum sample: 1. Use standard laboratory techniques for collection of venous (heparinized) blood samples. 2. Spin down using standard laboratory technique to separate plasma from blood cells if desired. Then open the sample blocking dilutor vial (lilac lid) and apply the correct volume of blood (20 µl), plasma or serum (8 µl), e.g. with a laboratory pipette. Close the sample blocking dilutor vial and mix gently by hand or vortex briefly.

Adding the sample: 1. Gloves should be used while performing the test. The test must be performed on a flat surface. 2. Open the sample blocking dilutor vial (lilac lid) containing the sample and completely withdraw its content with the syringe. Make sure no large air bubbles are formed inside. Place the syringe firmly on the fluid port of the assay device and inject the entire content of the syringe by pressing down the plunger. **Incubate for 8 minutes. Note!** If little backflow of the blood sample mix occurs, this does not influence the test results as long as the membrane is soaked with the blood sample mix.

Primary antibody incubation: 1. Open the first wash solution vial (blue lid), withdraw the solution completely and inject as described above. 2. Open the primary antibody vial (white lid, translucent tube), withdraw the solution completely and inject as described above. **Incubate for 8 minutes.**

Secondary antibody incubation: 1. Open the second wash solution vial (blue lid), withdraw the solution completely and inject as described above. 2. Open the secondary antibody vial (yellow lid), withdraw the solution completely and inject as described above. **Incubate for 8 minutes.**

Color reaction: 1. Open the third wash solution vial (blue lid), withdraw the solution completely and inject as described above. 2. Open the two color reagent vials (brown tubes, white lid), withdraw the solution of the first vial completely and inject as described above. Start the timer. **Within 30 sec:** Withdraw the solution of the second vial completely and inject as described above. **Incubate for 5 minutes** (total incubation time of the two steps).

Stopping the reaction: Open the stop buffer vial (green lid), withdraw the solution completely and inject as described above.

TEST PROCEDURE FOR TWO OR MORE PARALLEL TESTS

Experienced users of the AESKUCARE® Gluten Related Disorders IgA tests may perform more than one test in parallel. It is recommended to start the tests in intervals of at least 1 minute or multiples thereof to ensure proper handling during the development steps. It is not recommended to perform more than five tests in parallel.

INTERNAL CONTROLS

Each test row contains methodological controls (negative (NC) and function (FC) controls). FC signals appear as coloured bands after the test has been performed. The first and fourth test rows contain standards or cut-off calibrator that are essential for the quantitative and qualitative evaluation of the test results. Standard 1 in the first test row corresponds to the cut-off calibrator required for the qualitative evaluation of the test. **Test rows 1 and 4 (quantitative evaluation):** if standards or FC are missing or NC is \geq standard 1, the result for the whole test device will be invalid and the test will need to be repeated using a new test kit. **Test row 1 (qualitative evaluation):** if the cut-off calibrator or FC are missing or NC is \geq cut-off-calibrator, the result for the whole test device is invalid. **Test rows 2 and 3 (quantitative and qualitative evaluation):** if FC is missing or NC is \geq standard 1 / \geq cut-off calibrator, the results for the corresponding test row (test row 2 or 3) will be invalid and the test will need to be repeated using a new test kit.

ASSAY RESULT READING

Qualitative evaluation: To evaluate the results, position the test device with fluid port (sample inlet) pointing to the left. In this orientation the names of the antigens are positioned under the corresponding test areas. Results must be read in good light. Be cautious with shadows caused by the front label and avoid confusion with test signals. Start from the top left of the first test row (total IgA) and then read out from the left until the right end of second and third rows. This corresponds to the antigen layout in the evaluation sheet. For a qualitative evaluation of the results, compare the intensity of the cut-off calibrator (standard 1) in the first test row (shown with arrow on the upper label of the tests device) against the intensities of the individual antigens (visual average of triplicate measurement). If there is a significant difference between the signal intensities of an antigen within one triplicate, two signals with the most similar signal intensities should be selected for visual evaluation ( - ignore the weak signal;  - ignore the strong signal). Results are recorded in the corresponding evaluation sheet.

Possible assay results (qualitative evaluation)

negative	Signal intensity less than the cut-off calibrator (standard 1)
equivocal	Signal intensity approximately the same as the cut-off calibrator (standard 1)
positive	Signal intensity greater than the cut-off calibrator (standard 1)

Quantitative evaluation: If using a separately available reader, the results can be analysed quantitatively. To do so, place the test device in the reader and read out the results using the accompanying software and corresponding control chart. Note the instructions in the corresponding manual. The quantitative evaluation is performed automatically using a regression analysis in which the raw values of five standards (standards 0-4) are fitted against the concentration in U/ml (arbitrary units). A 4-parameter logistic function integrated into the control chart is used for the test analysis. The results for each antigen are expressed in U/ml.

Possible assay results (quantitative evaluation)

Result *	sIgA concentration (arbitrary units, U/ml)
negative	< 12.00
equivocal	12.00 – 17.99
low positive	18.00 – 29.99
positive	30.00 – 199.99
high positive	\geq 200.00

* refer to the section „Interpretation of Results“

Note 1. Take care to use only the evaluation sheet that is part of the corresponding AESKUCARE® Gluten Related Disorders IgA kit. 2. The reading of the test results should be done immediately after the test procedure to avoid confusion with other tests. Test results remain stable for at least 24 h if stored in the dark. The background color of the test areas might change with time. In case small air bubbles are present they may be removed by carefully inflating the device with some air using the syringe provided until the test areas are bubble free. 3. If the standard bands or cut-off calibrator are not visible, test results are invalid.

Analytical sensitivity has been assessed by multiple analysis of blank and low positive samples followed by calculation of the Limit of Blank (LoB) and Limit of Detection (LoD), respectively. The following overall values have been determined for all antigens: LoB = 0.24 U/ml; LoD = 0.87 U/ml.

Analytical specificity: AESKUCARE® Gluten Related Disorders IgA detects specific IgA antibodies. No significant influences on the test results could be observed with the interfering substances: triglycerides (2000 mg/dL), hemolysate (500 mg/dL), bilirubin (40 mg/dL), and IgG (4500 mg/dL). No cross-reactivity with other immunoglobulin species is known.

Limitations of the examination procedure: Clinically relevant interfering substances or limitations are not known. Carry-over effects are not known. Unsuitable primary samples: coagulated blood. Affecting factors/circumstances and precautions: refer to the sections "Sample Material", "Test Procedure", and "Assay Result Reading".

Relative sensitivity/specificity: The test results of AESKUCARE® Gluten Related Disorders IgA have been evaluated in comparison to the CE marked AESKUBLOTS® Gluten Related Disorders IgA test system (AESKU.DIAGNOSTICS GmbH & Co. KG). The comparison showed following values of relative sensitivity/specificity in percent: mTG neo 100/90, TG3 100/100, tTG 67/100, tTG neo 100/100, mTG -/100, Frazer's Fraction 75/95, gliadin 73/100, DGP 93/93.

* Previous studies have shown that the blood samples containing antibodies against mTG are quite rare. This observation has also been confirmed when using AESKUCARE® Gluten Related Disorders IgA. For this reason, the relative sensitivity of this antigen cannot be properly estimated.

Reference ranges: Due to the lack of an international reference standards, this assay is calibrated in arbitrary units (U/ml). For the threshold values used for positive result ranges of sIgA, refer to the section "Assay Result Reading".

Inter-assay precision: Three assays with control serum were performed by five persons to determine inter-assay precision (n = 15). The mean values were determined separately for each antigen.* The results are summarized in the table below. A decimal level difference of ≤ 1.0 between the assays is accepted for the inter-assay precision of all antigens.

* In statistical data analysis, the data obtained in U/ml have been categorized into decimal levels 0 – 3.5. In general, decimal levels ≤ 0.99 represent < 12.00 U/ml results, decimal levels ≥ 1.00 and ≤ 1.99 represent 12.00 – 29.99 U/ml results, decimal levels ≥ 2.00 and ≤ 2.99 represent 30.00 – 199.99 U/ml results, and decimal levels ≥ 3.00 and ≤ 3.50 represent 200.00 – 300.00 U/ml results.

Antigen	Number of measurements	Mean (decimal level)	Inter assay decimal level difference (max vs. min)	Analysis
total IgA	15	3.2	0.4	positive
mTG neo	15	2.6	0.2	positive
TG3	15	0.0	0.1	negative
tTG	15	2.5	0.4	positive
tTG neo	15	2.8	0.3	positive
mTG	15	0.0	0.0	negative
Frazer's F	15	2.8	0.4	positive
gliadin	15	2.4	0.2	positive
DGP	15	2.0	0.2	positive

Lot to Lot precision: Assay with control sera were performed with three lots. The Lot to Lot precision was 96%

Linearity: For chosen control sera, a linear correlation between dilution and antibody concentration has been shown. Correlation coefficient R² of the linear straight line: > 0.95 (90% of all positive antigens), > 0.94 (100% of all positive antigens). Due to the heterogeneity of human antibodies, it cannot be excluded that individual sera might show a non-linear behaviour.

QUALITY CONTROL

It is good laboratory practice to record lot numbers of all components used.

INTERPRETATION OF RESULTS

In the diagnostics of GRD it is not usual to consider particular values of sIgA as strictly defined limits. In general, with an increasing sIgA concentration toward GRD antigen-biomarkers, a direct correlation with symptoms is more probable. Signal intensities ≥ 30.0 U/ml correspond to a moderate or high sIgA level. **Note!** In order to make a nutritional recommendation or support a diagnosis of GRD, the test results must always be evaluated together with the medical history and any symptoms that occur.

GRDs are disorders caused by gluten, which include celiac disease (CD), non-celiac gluten sensitivity (NCGS), dermatitis herpetiformis (DH) and wheat allergy.

CD is a disorder of the small intestine with a chronic course. In affected individuals, the intake of gluten with food leads to a T-cell mediated immunological reaction.

Antibodies against gliadin, tTG and DGP correlate with the activity of the disease and thus are useful for professionals to support a diagnosis. Anti-gliadin antibodies may also be found in a group of healthy individuals.

Recent studies have shown that antibodies directed against gliadin of CD patients bind a very limited number of specific epitopes on the gliadin molecule. Selective deamidation of gliadin by tissue transglutaminase results in enhanced binding by anti-gliadin antibodies.

Assays using deamidated (DGP) and defined peptides have been shown to have higher diagnostic accuracy for CD when compared to standard anti-gliadin assays.

mTG shows substantial functional similarity with tTG in the process of selective deamidation of proteins. The cross-link of tTG or mTG with gliadin-specific peptides results in neo-epitopes of tTG and mTG, respectively. As these neo-epitopes are structurally closer to the physiological antigens, they show a markedly increased sensitivity and specificity. Furthermore, mTG neo-epitopes biomarkers may provide evidence not only for the early stages of CD but also for NCGS.

Frazer's fraction (also known as PT-gliadin) represents a peptic-tryptic digests of the gliadin, which retains the reactive properties of gliadin and mimic the fraction that reaches the duodenum after gastric digestion. The assumption is that a naturally occurring mixture of gliadin peptides significantly increases the sensitivity and specificity of anti-gliadin antibodies.

DH is a chronic, blistering skin disease associated with severe itching. DH patients possess antibodies against TG3. Additionally, granular IgA deposits can be detected under the skin. For all patients with DH, CD can usually be diagnosed as a comorbid condition. Therefore, the serological detection of the antibodies against specific biomarkers of CD (gliadin, Frazer's fraction, tTG, neo-epitopes of tTG and DGP) in combination with analysis of the sensitivity towards TG3 is a useful tool to support the diagnosis of DH.

The "total IgA" test parameter is used for the detection of possible IgA deficiency. Approximately 2% - 5% of CD patients display an IgA deficiency and thus negative IgA autoantibody tests (signal intensity of total IgA < 18 U/ml when the AESKUCARE® Gluten Related Disorders IgA test is evaluated quantitatively or lower than the cut-off calibrator when the test is evaluated qualitatively), even when the disease is active. The determination of specific IgG antibodies against GRD antigen-biomarkers is recommended for patients with confirmed IgA deficiency.

DISPOSAL

Precautions for the handling of human blood and safe waste disposal of used parts (containing blood) must be observed. Please refer to all local/national guidelines for biohazard handling. The rest of the kit may not be considered biohazard and can be disposed accordingly.

FURTHER LIMITATIONS

1. A final clinical diagnosis of GRD should not be based on the results of any single diagnostic method, but must only be made after all clinical and laboratory findings have been evaluated. In vitro evidence of sIgA against GRD antigen-biomarkers should always be accompanied by full medical history and analysis of any symptoms. 2. Negative in vitro results may occasionally occur in patients with GRD symptoms that clearly correlate with gluten-containing diet. 3. Patients with clinical symptoms suggestive of GRD but a negative result or low sIgA concentration should be referred to a specialist for further investigation. 4. Samples from apparently healthy individuals may also contain autoantibodies specific for GRDs. 5. It should be noted that under a gluten-free diet, specific antibodies against GRD antigen-biomarkers decrease significantly until they are no longer detectable. 6. Reliable and reproducible results will be obtained when the assay is performed according to the procedural instructions and when following good laboratory practices.

EXPECTED VALUES

To investigate the expected range, randomly selected donors have been tested on each of the GRD antigen-biomarkers.

Antigen	Number of measurements	Median [U/ml]	Mean [U/ml]	Standard Deviation [U/ml]
mTG neo	60	0.00	0.38	1.10
TG3	60	0.33	3.22	9.75
tTG	60	0.00	0.60	1.17
tTG neo	60	0.00	0.26	0.84
mTG	60	0.00	1.23	4.87
Frazer's F	60	0.00	0.48	1.12
gliadin	60	0.00	0.41	1.08
DGP	60	0.04	1.22	4.07

It is recommended that each laboratory establishes its own range of normal values.

In general, the expected values for apparently healthy individuals (children, men, and women) are < 18 U/ml at room temperature (18-30°C) for all GRD antigens of the test. For the threshold values of the positive result ranges of sIgA, refer to the section "Assay Result Reading".

WARRANTY

The performance data presented here was obtained using the procedure indicated. Any change or modification to the procedure not recommended here may affect the results, in which case AESKU. DIAGNOSTICS GmbH & Co. KG disclaims all warranties expressed, implied or statutory, including the implied warranty of merchantability and fitness for use. AESKU.DIAGNOSTICS GmbH & Co. KG and its authorized distributors, in such event, shall not be liable for damages indirect or consequential.

REFERENCES

1. Elli L, et al. Diagnosis of gluten related disorders: Celiac disease, wheat allergy and non-celiac gluten sensitivity. World J Gastroenterol. (2015) 21(23):7110-7119. 2. Caio G, et al. Celiac disease: a comprehensive current review. BMC Med. (2019) 17(1):142. 3. Frazer AC, et al. Gluten-induced enteropathy: the effect of partially digested gluten. Lancet (1959) 2(7097):252-255.

4. Lerner A, et al. Antibodies against neo-epitope tTG complexed to gliadin are different and more reliable than anti-tTG for the diagnosis of pediatric celiac disease. J Immunol Methods (2016) 429:15-20. 5. Lerner A, et al. Transglutaminases in dysbiosis as potential environmental drivers of autoimmunity. Front Microbiol. (2017) 8:66. 6. Lerner A and Matthias T. Possible association between celiac disease and bacterial transglutaminase in food processing: a hypothesis. Nutr Rev. (2015) 73(8):544-552. 7. Lytton SD, et al. Neo-epitope tissue transglutaminase autoantibodies as a biomarker of the gluten sensitive skin disease-dermatitis herpetiformis. Clin Chim Acta. (2013) 415:346-9. 8. Matthias T, et al. The industrial food additive microbial transglutaminase, mimics the tissue transglutaminase and is immunogenic in celiac disease patients. Autoimmun Rev. (2016) 15(12):1111-1119. 9. Mothes T. Deamidated gliadin peptides as targets for celiac disease-specific antibodies. Adv Clin Chem. (2007) 44:35-63. 10. Prause C, et al. Antibodies against deamidated gliadin as new and accurate biomarkers of childhood coeliac disease. J Pediatr Gastroenterol Nutr. (2009) 49(1):52-58. 11. Reunala T, et al. Dermatitis herpetiformis: A common extraintestinal manifestation of coeliac disease. Nutrients (2018) 10(5): 602. 12. Salmi TT. Dermatitis herpetiformis. Clin Exp Dermatol. (2019) 44(7):728-731. 13. Singh P. Global prevalence of celiac disease: systematic review and meta-analysis. Clin Gastroenterol Hepatol. (2018) Jun;16(6):823-836.

Guidelines: CLSI GP44-A4: Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests. CLSI EP17-A2: Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures.

PATENTS/TRADEMARKS




EU-Patent: EP 2 839 290 B1, US-Patent: US 10,571,466 B2

ADDRESS

Manufactured by AESKU.DIAGNOSTICS GmbH & Co. KG
Mikroforum Ring 2, 55234 Wendelsheim, Germany
Tel: +49 6734 - 9622-0, Fax: +49 6734 - 9622-2222
Web: www.aesku.com

Additional company and product information can be requested at any time.

The instructions for use apply to AESKUCARE® Gluten Related Disorders IgA (REF: 860001).

	In vitro diagnostic medical device
	Consult instructions for use
	Article number
	Batch code
	Use by
	Content sufficient for < 1 > test
	Do not reuse
	Temperature limitation
	Caution
	CE Marking of conformity
	Manufacturer

HISTORY: Removed information of external quality controls. Section of internal controls extended.

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REFERENZEN

1. Elli L, et al. Diagnosis of gluten related disorders: Celiac disease, wheat allergy and non-celiac gluten sensitivity. World J Gastroenterol. (2015) 21(23):7110-7119.
2. Caio G, et al. Celiac disease: a comprehensive current review. BMC Med. (2019) 17(1):142.
3. Frazer AC, et al. Gluten-induced enteropathy: the effect of partially digested gluten. Lancet (1959) 2(7097):252-255.
4. Lerner A, et al. Antibodies against neo-epitope tTG complexed to gliadin are different and more reliable than anti-tTG for the diagnosis of pediatric celiac disease. J Immunol Methods (2016) 429:15-20.
5. Lerner A, et al. Transglutaminases in dysbiosis as potential environmental drivers of autoimmunity. Front Microbiol. (2017) 8:66.
6. Lerner A and Matthias T. Possible association between celiac disease and bacterial transglutaminase in food processing: a hypothesis. Nutr Rev. (2015) 73(8):544-552.
7. Lytton SD, et al. Neo-epitope tissue transglutaminase autoantibodies as a biomarker of the gluten sensitive skin disease-dermatitis herpetiformis. Clin Chim Acta. (2013) 415:346-9.
8. Matthias T, et al. The industrial food additive microbial transglutaminase, mimics the tissue transglutaminase and is immunogenic in celiac disease patients. Autoimmun Rev. (2016) 15(12):1111-1119.
9. Mothes T. Deamidated gliadin peptides as targets for celiac disease-specific antibodies. Adv Clin Chem. (2007) 44:35-63.
10. Prause C, et al. Antibodies against deamidated gliadin as new and accurate biomarkers of childhood coeliac disease. J Pediatr Gastroenterol Nutr. (2009) 49(1):52-58.
11. Reunala T, et al. Dermatitis herpetiformis: A common extraintestinal manifestation of coeliac disease. Nutrients (2018) 10(5): 602.
12. Salmi TT. Dermatitis herpetiformis. Clin Exp Dermatol. (2019) 44(7):728-731.
13. Singh P. Global prevalence of celiac disease: systematic review and meta-analysis. Clin Gastroenterol Hepatol. (2018) Jun;16(6):823-836.

Richtlinien: CLSI GP44-A4: Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests. CLSI EP17-A2: Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures.

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

REFERENCIAS

1. Elli L, et al. Diagnosis of gluten related disorders: Celiac disease, wheat allergy and non-celiac gluten sensitivity. World J Gastroenterol. (2015) 21(23):7110-7119. 2. Cao G, et al. Celiac disease: a comprehensive current review. BMC Med. (2019) 17(1):142. 3. Frazer AC, et al. Gluten-induced enteropathy: the effect of partially digested gluten. Lancet (1959) 2(7097):252-255. 4. Lerner A, et al. Antibodies against neo-epitope tTG complexed to gliadin are different and more reliable than anti-tTG for the diagnosis of pediatric celiac disease. J Immunol Methods (2016) 429:15-20. 5. Lerner A, et al. Transglutaminases in dysbiosis as potential environmental drivers of autoimmunity. Front Microbiol. (2017) 8:666. 6. Lerner A and Matthias T. Possible association between celiac disease and bacterial transglutaminase in food processing: a hypothesis. Nutr Rev. (2015) 73(8):544-552. 7. Lytton SD, et al. Neo-epitope tissue transglutaminase autoantibodies as a biomarker of the gluten sensitive skin disease-dermatitis herpetiformis. Clin Chim Acta. (2013) 415:346-9. 8. Matthias T, et al. The industrial food additive microbial transglutaminase, mimics the tissue transglutaminase and is immunogenic in celiac disease patients. Autoimmun Rev. (2016) 15(12):1111-1119. 9. Mothes T. Deamidated gliadin peptides as targets for celiac disease-specific antibodies. Adv Clin Chem. (2007) 44:35-63. 10. Prause C, et al. Antibodies against deamidated gliadin as new and accurate biomarkers of childhood coeliac disease. J Pediatr Gastroenterol Nutr. (2009) 49(1):52-58. 11. Reunala T, et al. Dermatitis herpetiformis: A common extraintestinal manifestation of coeliac disease. Nutrients (2018) 10(5): 602. 12. Salmi TT. Dermatitis herpetiformis. Clin Exp Dermatol. (2019) 44(7):728-731. 13. Singh P, Global prevalence of celiac disease: systematic review and meta-analysis. Clin Gastroenterol Hepatol. (2018) Jun;16(6):823-836.

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4. Lerner A, et al. Antibodies against neo-epitope tTG complexed to gliadin are different and more reliable than anti-tTG for the diagnosis of pediatric celiac disease. J Immunol Methods (2016) 429:15-20. 5. Lerner A, et al. Transglutaminases in dysbiosis as potential environmental drivers of autoimmunity. Front Microbiol. (2017) 8:66. 6. Lerner A and Matthias T. Possible association between celiac disease and bacterial transglutaminase in food processing: a hypothesis. Nutr Rev. (2015) 73(8):544-552. 7. Lytton SD, et al. Neo-epitope tissue transglutaminase autoantibodies as a biomarker of the gluten sensitive skin disease-dermatitis herpetiformis. Clin Chim Acta. (2013) 415:346-9. 8. Matthias T, et al. The industrial food additive microbial transglutaminase, mimics the tissue transglutaminase and is immunogenic in celiac disease patients. Autoimmun Rev. (2016) 15(12):1111-1119. 9. Mothes T. Deamidated gliadin peptides as targets for celiac disease-specific antibodies. Adv Clin Chem. (2007) 44:35-63. 10. Prause C, et al. Antibodies against deamidated gliadin as new and accurate biomarkers of childhood coeliac disease. J Pediatr Gastroenterol Nutr. (2009) 49(1):52-58. 11. Reunala T, et al. Dermatitis herpetiformis: A common extraintestinal manifestation of coeliac disease. Nutrients (2018) 10(5): 602. 12. Salmi TT. Dermatitis herpetiformis. Clin Exp Dermatol. (2019) 44(7):728-731. 13. Singh P, Global prevalence of celiac disease: systematic review and meta-analysis. Clin Gastroenterol Hepatol. (2018) Jun; 16(6):823-836.

Linee guida: CLSI GP44-A4: Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests. CLSI EP17-A2: Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures.

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Prodotto da AESKU.DIAGNOSTICS GmbH & Co. KG