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REFERENZEN: 1. Williams PB, et al. J. Efficacy of a single diagnostic test for sensitisation to common inhalant allergens. Ann Allergy Asthma Immunol 2001;86:196-202. 2. Ahlstedt S. Understanding the usefulness of specific IgE blood tests in allergy. Clin Exp Allergy 2002;32:11-6. 3. Martinez FD, et al. Asthma and wheezing in the first six years of life. N Engl J Med 1995;332:133-8. 4. Gern JE, Busse WW. The role of viral infections in the natural history of asthma. J Allergy Clin Immunol 2000;106:201-12. 5. Heinzerling L et al. Standard skin prick testing and sensitization to inhalant allergens across Europe – a survey from the GA2LEN network. Allergy 2005;60:1287-1300. 6. Wever-Hess J, et al., Prognostic characteristics of asthma diagnosis in early childhood in clinical practice. Acta Paediatr 1999;88:827-34. 7. Eysink PED, et al. Development of specific immunoglobulin E in coughing toddlers: A medical records review of symptoms in general practice. Pediatr Allergy Immunol 2001;12:133-141. 8. Duran-Tauleria E, et al. The utility of specific immunoglobulin E measurements in primary care. Allergy 2004;59 (Suppl.78):35-41.

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INSTRUCTIONS FOR USE (EN)

INTENDED USE: Aeskucare® Allergy is an in vitro semi-quantitative enzyme immunoassay for the parallel measurement of allergen-specific human IgE antibodies to 20 allergens and allergen mixtures in heparinized or Na-EDTA blood, venous or capillary blood, plasma or serum. The test is intended to be a point of care tool for comprehensive first line in vitro allergy diagnostics in conjunction with other clinical findings, and is to be carried out by healthcare professionals experienced in the use of in vitro diagnostic methods.

PRINCIPLE OF THE PROCEDURE: Aeskucare® Allergy is a modified ELISA test system (Enzyme Linked Immunosorbent Assay). The assay device contains 4 rows with 5 test areas each. Each test area comprises a specific allergen or allergen mixture and two controls. Specific IgE antibodies (sIgE) from the sample bind to the corresponding antigens of the allergen. The bound human IgE is detected by specific anti-human-IgE antibodies. A secondary antibody coupled to alkaline phosphatase is added for signal amplification. After adding color reagent, any bound IgE is identified by an indigo-colored precipitate. For all 20 allergens, individual sensitization levels can be read semi-quantitatively in 5 levels with the naked eye. The read-out is performed by comparing the color saturation (intensity) of the signal with those of the corresponding controls in each test area.

Precautions: 1. For human in vitro diagnostic use only. 2. Please read the entire contents of these Instructions for Use before performing the test. 3. Do not use reagents beyond their expiration dates. 4. In case of damage to the packaging, please inspect the protecting bag of the test device carefully for damage and make sure that the protection of the fluid port is intact. Make sure the reagents are not damaged or open. In case of doubt do not use the test kit in order to avoid incorrect results or a wrong diagnosis. 5. Do not pool any reagents. 6. The color reagent contains 5-bromo-4-chloro-3-indolylphosphate and nitro-blue tetrazolium chloride. Avoid contact with skin. Wear suitable protection gloves. For more information see Material Safety Data Sheet (available on request). 7. All components are for single use only. 8. The test results are invalid, if the precautions are not followed. 9. Fluid paths can be blocked by coagulated blood.

STORAGE: Store in a dark and cool place at 2-8°C/36-46°F. All reagents are ready-to-use and packaged in screw tubes. **EXPIRATION:** The expiration date is given on the outer packaging. The expiration date of the kit is valid for all kit components, even if the expiration of single components is different! After expiration, all test components have to be discarded.

REAGENTS AND MATERIAL: 1. **DILUT - Sample dilutor**, 300 µL, 1x (red lid) 2. **WASH - Washing solution**, 1 mL each, 3x (blue lid) 3. **Ab1 - Anti body solution 1**, 800 µL, 1x (grey lid) 4. **Ab2 - Anti body solution 2**, 800 µL, 1x (yellow lid) 5. **COLR - Color reagent** 800 µL each (contains 5-bromo-4-chloro-3-indolylphosphate/nitro-blue tetrazolium chloride), 2x (white lid) 6. **STOP - Stop buffer**, 1 mL, 1x (green lid). **Note!** All reagents are ready to use; Store at 2-8 °C until the expiration date.

Human material used in the Aeskucare® Allergy assay device components was tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as being potentially infectious.

Materials provided: instructions for use, evaluation sheet Aeskucare® Allergy assay device, reagent rack containing 9 color-coded vials, blood sampling pouch: alcohol wipe, CE, safety lancet sterile CE 0124, 200 µL capillary CE with lithium heparin, gauze / absorbent pad CE, band-aid CE, syringe **Additional material, not provided in the kit:** disposable gloves, timer. For venous blood and/or plasma samples: normal venipuncture supplies, heparinized blood collection tubes, laboratory pipette. **List of allergens and allergen extracts:** see evaluation sheet in the packaging.



TREATMENT OF TEST SAMPLES: Capillary blood from a finger prick, from earlobes or heels (suitable for infants) can be used in this test. The blood sampling capillary provided in the kit should be used to collect the correct volume (200 µL). Heparinized (sodium or lithium heparin) venous whole blood (200 µL), plasma or serum (67 µL) can also be used as samples in this test.

Collection of venous blood, plasma, or serum sample: 1. Use normal laboratory techniques for collecting venous (heparinized) blood samples. 2. If desired, centrifuge the blood according to normal laboratory technique to separate plasma from blood cells. Laboratory pipettes should be used to dispense the required volume of plasma or serum. Close the sample dilutor vial and mix carefully by hand or by shortly vortexing the sample. Once the blood sample has been added to the sample dilutor, all samples may be stored at 2-8 °C for up to 72 hours. Make sure to label the test kit with the patient ID, if you want to store the sample. **Note!** Do not store the sample in the blood sampling capillary.

TEST PROCEDURE

Preparation: All reagents and the assay device must be adjusted to room temperature (18-30°C) for at least 30 minutes before use. 1. Open the kit and take out the assay device, the syringe, the reagent rack, and the blood sampling pouch. Place them on a flat surface and 2. Open and assemble the reagent rack. 3. Label the assay device with the patient ID in the grey labelling area.

Collection of capillary blood sample: Open the blood sampling pouch and take out its contents: alcohol wipe, safety lancet, 200 µL blood sampling capillary, gauze, and band aid. Place all items on a clean, flat surface. 1. Massage your fingertip or put it into warm water for a few moments before collecting the blood so as to optimize the blood flow. 2. Clean the site to be pricked with the alcohol wipe. 3. Unlock the safety lancet as shown (fig. step A) and pull off the protective cap. 4. Press the safety lancet to the finger site to be pricked (fig. step B) and firmly press the button.

Note! Make sure the hand is below the heart level. If the blood flow is initially slow, gently massage the finger and check if the hand is below (!) the heart.

1. Collect the drops of blood with the tip of the blood sampling capillary. Hold the blood sampling capillary horizontally or at a slightly downward angle and let the blood flow until the capillary is properly filled. Make sure not to cover the ventilation hole in the plunger. The filling mark shows the sample volume (200 µL). 2. Cover the wound with the band-aid. 3. Immediately afterwards open the sample dilutor vial (red lid) and put the tip of the blood sampling capillary in a position just below the vial opening. Cover the hole at the end of the plunger with your finger. Gently press the plunger to transfer the blood into the sample dilutor vial. Close the sample dilutor vial and mix gently by hand or vortex briefly.

Note! Do not store the blood without having added it to the sample dilutor. Once added to the sample dilutor, all samples may be stored at 2-8 °C for up to 72 hours. If untreated blood is used (i.e. without adding the sample dilutor) then the tests results are not valid.

Adding the sample: 1. Gloves should be used while performing the test. 2. Remove the sticker from the fluid port of the assay device. 3. Open the sample dilutor vial (red lid) containing the sample and completely withdraw its content with the syringe. Make sure no large air bubbles are formed inside. Place the syringe firmly on the fluid port of the assay device and inject the entire content of the syringe by pressing down the plunger. **Incubate for 8 minutes.** Note! If little backflow of the blood sample mix occurs, this does not influence the test results as long as the membrane is soaked with the blood sample mix.

Primary antibody incubation: 1. Open the first vial with washing solution (WASH - blue lid), completely withdraw the solution and inject it as described above. 2. Open the primary antibody vial (Ab1 - transparent lid), completely withdraw the solution and inject it as described above. **Incubate for 8 minutes.**

Secondary antibody incubation: 1. Open a second washing solution vial (WASH - blue lid), completely withdraw the solution and inject it as described above. 2. Open the secondary antibody vial (Ab2 - yellow lid), completely withdraw the solution and inject it as described above. **Incubate for 8 minutes.**

Colour reagent: 1. Open a third washing solution vial (WASH - blue lid), completely withdraw the solution and inject it as described above. 2. Open the two color reagent vials (COLR - brown tubes, white lid), completely withdraw the solution of the first vial and inject it as described above. Start the timer. 3. **Within 30 sec**, completely withdraw the solution of the second vial and inject it as described above. **Incubate for 5 minutes** (total incubation time of the two steps).

Stopping the reaction: 1. Open the stop buffer vial (STOP - green lid), completely withdraw the solution and inject it as described above. 2. Cover the fluid port with the sticker (optional).

TEST PROCEDURE FOR TWO OR MORE PARALLEL TESTS: Experienced users of Aeskucare® Allergy tests may perform several tests in parallel. It is recommended to test the tests in intervals of at least 1 min or multiples thereof to ensure proper handling during the development steps. It is not recommended to perform more than five tests in parallel.

ASSAY RESULT READING: To evaluate the results, position the assay device with the fluid port (sample inlet) pointing to the left. In this orientation allergen codes are positioned under the corresponding test area. Results must be read in good light conditions. Make sure to prevent shadowing. Start from the top left and then read each test area until the right end of the row. Repeat row by row. This corresponds to the allergen layout in the evaluation sheet. To estimate the level (1-5) of a signal band, its intensity must be compared with the intensity of the two procedural control bands (lower and higher standard) on the left and right sides of the corresponding test area. All 20 results are recorded on the corresponding evaluation sheet.

Note! 1. Take care to use only evaluation sheets taken from the corresponding Aeskucare® Allergy kit. 2. The reading of the test results should be done immediately after the test procedure to avoid confusion with other tests. Test results remain stable for at least 24 h if stored in the dark. The background of the results in the test fields might change over time. In case small air bubbles have developed in the test device, they may be removed by carefully inflating the device with some air using the syringe if the control bands are not visible, the test results are invalid.

Possible assay evaluation results

Level 1	Signal intensity less than the lower standard	
Level 2	Signal intensity approximately the same as lower standard	
Level 3	Signal intensity greater than the lower and less than the higher standard	
Level 4	Signal intensity approximately the same as higher standard	
Level 5	Signal intensity greater than the higher standard	

Correlation between Aeskucare® Allergy-Level, sensitization, class and sIgE*			
Aeskucare® Allergy-Level	Sensitization	CAP Class	sIgE [kU _e /L]
1	low	0-2	0-3.4
2*	moderate	2-3	0.7-17.4
3	high	3-4	3.5-49.9
4*	very high	4-5	17.5-100
5	very high	5-6	> 50

* The two Aeskucare® Allergy controls are identical to levels 2 and 4. CAP class 2 is usually seen as a borderline sensitization.

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REFERENCIAS: 1. Williams PB, et al. J. Efficacy of a single diagnostic test for sensitisation to common inhalant allergens. Ann Allergy Asthma Immunol 2001;86:196-202. 2. Ahlstedt S. Understanding the usefulness of specific IgE blood tests in allergy. Clin Exp Allergy 2002;32:11-6. 3. Martinez FD, et al. Asthma and wheezing in the first six years of life. N Engl J Med 1995;332:133-8. 4. Gern JE, Busse WW. The role of viral infections in the natural history of asthma. J Allergy Clin Immunol 2000;106:201-12. 5. Heinzlering L et al. Standard skin prick testing and sensitization to inhalant allergens across Europe – a survey from the GAZLEN network. Allergy 2005;60:1287-1300. 6. Wever-Hess J, et al. Prognostic characteristics of asthma diagnosis in early childhood in clinical practice. Acta Paediatr 1999;88:827-34. 7. Eysink PED, et al. Development of specific immunoglobulin E in coughing toddlers: A medical records review of symptoms in general practice. Pediatr Allergy Immunol 2001;12:133-141. 8. Duran-Tauleira E, et al. The utility of specific immunoglobulin E measurements in primary care. Allergy 2004;59 (Suppl.78):35-41.

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
BIBLIOGRAFIA: 1. Williams PB, et al. J. Efficacy of a single diagnostic test for sensitisation to common inhalant allergens. Ann Allergy Asthma Immunol 2001;86:196-202. 2. Ahlstedt S. Understanding the usefulness of specific IgE blood tests in allergy. Clin Exp Allergy 2002;32:11-6. 3. Martinez FD, et al. Asthma and wheezing in the first six years of life. N Engl J Med 1995;332:133-8. 4. Gem JE, Busse WW. The role of viral infections in the natural history of asthma. J Allergy Clin Immunol 2000;106:201-12. 5. Heinzerling L et al. Standard skin prick testing and sensitization to inhalant allergens across Europe - a survey from the GA2LEN network. Allergy 2005;60:1287-1300. 6. Wever-Hess J, et al.. Prognostic characteristics of asthma diagnosis in early childhood in clinical practice. Acta Paediatr 1999;88:827-34. 7. Eysink PED, et al. Development of specific immunoglobulin E in coughing toddlers: A medical records review of symptoms in general practice. Pediatr Allergy Immunol 2001;12:133-141. 8. Duran-Tauleria E, et al. The utility of specific immunoglobulin E measurements in primary care. Allergy 2004;59 (Suppl.78):35-41.

INDIRIZZO: Realizzato da AESKU.Diagnostics GmbH & Co. KG, Mikroforum Ring 2, 55234 Wendelsheim, Germany, Tel: +49 6734 - 9622-0, Fax: +49 6734 - 9622-2222, Web: www.aesku.com

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ODKAZY: 1. Williams PB, et al. J. Efficacy of a single diagnostic test for sensitisation to common inhalant allergens. Ann Allergy Asthma Immunol 2001;86:196-202. 2. Ahlstedt S. Understanding the usefulness of specific IgE blood tests in allergy. Clin Exp Allergy 2002;32:11-6. 3. Martinez FD, et al. Asthma and wheezing in the first six years of life. N Engl J Med 1995;332:133-8. 4. Gern JE, Busse WW. The role of viral infections in the natural history of asthma. J Allergy Clin Immunol 2000;106:201-12. 5. Heinzerling L et al. Standard skin prick testing and sensitization to inhalant allergens across Europe – a survey from the GA2LEN network. Allergy 2005;60:1287-1300. 6. Wever-Hess J, et al.. Prognostic characteristics of asthma diagnosis in early childhood in clinical practice. Acta Paediatr 1999;88:827-34. 7. Eysink PED, et al. Development of specific immunoglobulin E in coughing toddlers: A medical records review of symptoms in general practice. Pediatr Allergy Immunol 2001;12:133-141. 8. Duran-Tauleria E, et al. The utility of specific immunoglobulin E measurements in primary care. Allergy 2004;59 (Suppl.78):35-41.

ADRESA VÝROBCU: AESKU.Diagnostics GmbH & Co. KG, Mikroforum Ring 2, 55234 Wendelsheim, Germany, Tel: +49 6734 - 9622-0, Fax: +49 6734 - 9622-2222, Web: www.aesku.com

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