AESKUCARE® IgA

Gluten Related Disorders Evaluation sheet (EN)

LOT XX.XX.XX
XXXX-XX

test row 1
IgA Std: 4 3 2 1 0 FC NC

test row 2
mTG neo TG3 tTG tTG neo FC NC

test row 3
mTG Frazer's F Gliadin DGP FC NC

test row 4
IgA Std: 4 3 2 1 0 FC NC

INVO

AESKUCARE®
Gluten Related Disorders
IgA

For visual evaluation of the test, compare the intensities of the individual antigens (visual average of three bands measurement) with the intensity of the cut-off calibrator shown by the red color. If there is a significant difference between the signal intensities of an antigen within one triplicate, two signals with the most similar signal intensities should be selected for visual evaluation.

Example: - ignore the weak signal; - ignore the strong signal

Test row 1: If the cut-off calibrator or the function control (FC) is missing or the negative control (NC) is \geq cut-off, the result for the whole test device is invalid. **Test rows 2 and 3:** If the function control (FC) is missing or the negative control (NC) is \geq cut-off, the results for the relevant test row are invalid.

Antigen	< cut-off	\geq cut-off	IgA deficiency: stop antigen evaluation. Determination of IgG antibodies against antigens specific for gluten related disorders is required. IgA normal level: proceed with the evaluation of the test.
total IgA (test row 1)	<input type="checkbox"/>	<input type="checkbox"/>	
Results	deficiency	normal level	

Antigen signal	Test result*
Signal intensity less than the cut-off	negative
Signal intensity approximately the same as the cut-off	equivocal
Signal intensity greater than the cut-off	positive

* These results are regarded clinical relevant if they are in correlation with the medical history and symptoms.

! **Make sure to read under proper light conditions!**

(Light grey test areas are not relevant for the visual evaluation of the test)

Patient: _____

Date of birth: _____

Date of test: _____

Signature: _____


Medical history: _____

Antigen	< cut-off	= cut-off	> cut-off
mTG neo microbial transglutaminase (neo-epitopes) (test row 2)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
TG3 transglutaminase 3 (human epidermal transglutaminase) (test row 2)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
tTG human tissue transglutaminase (test row 2)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
tTG neo human tissue transglutaminase (neo-epitopes) (test row 2)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
mTG microbial transglutaminase (test row 3)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Frazer's F Frazer's fraction (test row 3)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Gliadin alpha-gliadin (test row 3)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
DGP deamidated gliadin-specific peptide (test row 3)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Result	negative	equivocal	positive

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