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

## 11 Bibliografia

1. Y. Okada, H. Nagase, and E. D. Harris, Jr. A metalloproteinase from human rheumatoid synovial fibroblasts that digests connective tissue matrix components. Purification and characterization. *J.Biol.Chem.* 261 (30):14245-14255, 1986.
2. K. Obata, K. Iwata, Y. Okada, Y. Kohrin, E. Ohuchi, S. Yoshida, M. Shinmei, and T. Hayakawa. A one-step sandwich enzyme immunoassay for human matrix metalloproteinase 3 (stromelysin-1) using monoclonal antibodies. *Clin.Chim.Acta* 211 (1-2):59-72, 1992.
3. J. Martel-Pelletier, R. McCollum, N. Fujimoto, K. Obata, JM. Cloutier, JP. Pelletier. Excess of metalloproteases over tissue inhibitor of metalloprotease may contribute to cartilage degradation in osteoarthritis and rheumatoid arthritis. *Lab Invest.* 70 (6):807-15, 1994.
4. S. Sasaki, H. Iwata, N. Ishiguro, K. Obata, and T. Miura. Detection of stromelysin in synovial fluid and serum from patients with rheumatoid arthritis and osteoarthritis. *Clin.Rheumatol.* 13 (2):228-233, 1994.
5. Y. Yoshihara, K. Obata, N. Fujimoto, K. Yamashita, T. Hayakawa, and M. Shimmei. Increased levels of stromelysin-1 and tissue inhibitor of metalloproteinases-1 in sera from patients with rheumatoid arthritis. *Arthritis Rheum.* 38 (7):969-975, 1995.
6. H. Yamanaka, Y. Matsuda, M. Tanaka, W. Sendo, H. Nakajima, A. Taniguchi, and N. Kamatani. Serum matrix metalloproteinase 3 as a predictor of the degree of joint destruction during the six months after measurement, in patients with early rheumatoid arthritis. *Arthritis Rheum.* 43 (4):852-858, 2000.
7. M. D. Posthumus, P. C. Limburg, J. Westra, M. A. Van Leeuwen, and M. H. van Rijswijk. Serum matrix metalloproteinase 3 levels during treatment with sulfasalazine or combination of methotrexate and sulfasalazine in patients with early rheumatoid arthritis. *J.Rheumatol.* 29 (5):883-889, 2002.
8. Katrib, M. D. Smith, M. J. Ahern, J. Slavotinek, L. Stafford, C. Cuello, J. V. Bertouch, H. P. McNeil, and P. P. Youssef. Reduced chemokine and matrix metalloproteinase expression in patients with rheumatoid arthritis achieving remission. *J.Rheumatol.* 30 (1):10-21, 2003.
9. Tchvetverikov, L. R. Lard, J. DeGroot, N. Verzijl, J. M. TeKoppele, F. C. Breedveld, T. W. Huizinga, and R. Hanemaaijer. Matrix metalloproteinases-3, -8, -9 as markers of disease activity and joint damage progression in early rheumatoid arthritis. *Ann.Rheum.Dis.* 62 (11):1094-1099, 2003.
10. M. J. Green, A. K. Gough, J. Devlin, J. Smith, P. Astin, D. Taylor, and P. Emery. Serum MMP-3 and MMP-1 and progression of joint damage in early rheumatoid arthritis. *Rheumatology.(Oxford)* 42 (1):83-88, 2003.
11. Litinsky, D. Paran, D. Levartovsky, I. Wigler, I. Kaufman, I. Yaron, M. Yaron, D. Caspi, O. Elkayam. The effects of leflunomide on clinical parameters and serum levels of IL-6, IL-10, MMP-1 and MMP-3 in patients with resistant rheumatoid arthritis. *Cytokine.* 33 (2):106-10, 2006.
12. S. Visvanathan, M. U. Rahman, E. Keystone, M. Genovese, L. Klareskog, E. Hsia, M. Mack, J. Buchanan, M. Elashoff, and C. Wagner. Association of serum markers with improvement in clinical response measures after treatment with golimumab in patients with active rheumatoid arthritis despite receiving methotrexate: results from the GO-FORWARD study. *Arthritis Res.Ther.* 12 (6):R211, 2010.
13. Mahemara, T. Sugimoto, D. Sugiyama, S. Morinobu, G. Tsuji, S. Kawano, A. Morinobu, and S. Kumagai. Serum Matrix Metalloproteinase-3 as Predictor of Joint Destruction in Rheumatoid Arthritis, Treated with Non-biological Disease Modifying Anti-Rheumatic Drugs. *Kobe J. Med. Sci.* 56 (3):98-107, 2010.
14. Y. Urata, R. Uesato, D. Tanaka, Y. Nakamura, and S. Motomura. Treating to target matrix metalloproteinase 3 normalisation together with disease activity score below 2.6 yields better effects than each alone in rheumatoid arthritis patients: T-4 Study. *Ann.Rheum.Dis.*, 2011.
15. N. Nishimoto, K. Amano, Y. Hirabayashi, T. Horiuchi, T. Ishii, M. Iwahashi, M. Iwamoto, H. Kohsaka, M. Kondo, T. Matsubara, T. Mimura, H. Miyahara, S. Ohta, Y. Saeki, K. Saito, H. Sano, K. Takasugi, T. Takeuchi, S. Tohma, T. Tsuru, Y. Ueki, J. Yamana, J. Hashimoto, T. Matsutani, M. Murakami, and N. Takagi. Drug free REmission/low disease activity after cessation of tocilizumab (Actemra) Monotherapy (DREAM) study. *Mod.Rheumatol.*, 2013.
16. Hiura, S. Iwaki-Egawa, T. Kawashima, S. I. Fujisawa, T. Takeda, H. Komori, and Y. Watanabe. The diagnostic utility of matrix metalloproteinase-3 and high-sensitivity C-reactive protein for predicting rheumatoid arthritis in anti-cyclic citrullinated peptide antibody-negative patients with recent-onset undifferentiated arthritis. *Rheumatol.Int.*, 2013.

## 12 Simboli normativi

|                                                                                     |                                                                                                                                                                                                             |                                                                                                                                        |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| <b>IVD</b>                                                                          | - Diagnosi in vitro<br>- Pour diagnostic in vitro<br>- In Vitro Diagnostikum<br>- Para uso Diagnóstico in vitro                                                                                             | - For in vitro diagnostic use<br>- Para uso diagnóstico in vitro<br>- In Vitro Διαγνωστικό μέσο                                        |
| <b>REF</b>                                                                          | - Numero d'ordine<br>- Référence Catalogue<br>- Bestellnummer<br>- Número de catálogo                                                                                                                       | - Catalogue number<br>- Número de catálogo<br>- Αριθμός παραγγελίας                                                                    |
| <b>LOT</b>                                                                          | - Descrizione lotto<br>- Lot<br>- Chargen Bezeichnung<br>- Lote                                                                                                                                             | - Lot<br>- Lote<br>- Χαρακτηρισμός παρτίδας                                                                                            |
| <b>UDI</b>                                                                          | - Identificatore univoco del dispositivo<br>- Identifiant unique de l'appareil<br>- eindeutige Produktidentifizierung<br>- Identificador único do dispositivo                                               | - Unique Device Identifier<br>- Identificador único del dispositivo<br>- Μοναδικό αναγνωριστικό συσκευής                               |
| <b>CE</b>                                                                           | - Conformità europea<br>- Déclaration CE de Conformité<br>- Europäische Konformität<br>- Declaração CE de Conformidade                                                                                      | - EC Declaration of Conformity<br>- Declaración CE de Conformidad<br>- Ευρωπαϊκή συμφωνία                                              |
|    | - 96 determinazioni<br>- 96 tests<br>- 96 Bestimmungen<br>- 96 Testes                                                                                                                                       | - 96 tests<br>- 96 pruebas<br>- 96 προσδιορισμοί                                                                                       |
|    | - Rispettare le istruzioni elettroniche per l'uso<br>- Voir les instructions d'utilisation électronique<br>- Elektronische Gebrauchsanweisung beachten<br>- Seguir as instruções electrónicas de utilização | - See electronic instructions for use<br>- Siga las instrucciones electrónicas de uso<br>- Ακολουθήστε τις ηλεκτρονικές οδηγίες χρήσης |
|   | - Da utilizzarsi entro<br>- Utilise avant le<br>- Verwendbar bis<br>- Utilizar antes de                                                                                                                     | - Use by<br>- Utilizar antes de<br>- Χρήση μέχρι                                                                                       |
|  | - Conservare a 2-8°C (35,6-46,4°F)<br>- Conserver à 2-8°C (35,6-46,4°F)<br>- Lagerung bei 2-8°C (35,6-46,4°F)<br>- Conservar entre 2-8°C (35,6-46,4°F)                                                      | - Store at 2-8°C (35,6-46,4°F)<br>- Conservar a 2-8°C (35,6-46,4°F)<br>- Φυλάσσεται στους 2-8°C (35,6-46,4°F)                          |
|  | - Prodotto da<br>- Fabriqué par<br>- Hergestellt von<br>- Fabricado por                                                                                                                                     | - Manufactured by<br>- Fabricado por<br>- Κατασκευάζεται από                                                                           |
| <b>CO-CAL</b>                                                                       | - Calibratore cut-off<br>- Etalon Seuil<br>- Grenzwert Kalibrator<br>- Calibrador de cut-off                                                                                                                | - Cut off Calibrator<br>- Calibrador de cut-off<br>- Οριακός ορός Αντιδραστήριο βαθμονόμησης                                           |
| <b>CON +</b>                                                                        | - Controllo positivo<br>- Contrôle Positif<br>- Positiv Kontrolle<br>- Controllo positivo                                                                                                                   | - Positive Control<br>- Control Positivo<br>- Θετικός ορός ελέγχου                                                                     |
| <b>CON -</b>                                                                        | - Controllo negativo<br>- Contrôle Négatif<br>- Negativ Kontrolle<br>- Controllo negativo                                                                                                                   | - Negative Control<br>- Control Negativo<br>- Αρνητικός ορός ελέγχου                                                                   |
| <b>CAL</b>                                                                          | - Calibratore<br>- Etalon<br>- Kalibrator<br>- Calibrador                                                                                                                                                   | - Calibrator<br>- Calibrador<br>- Αντιδραστήριο βαθμονόμησης                                                                           |
| <b>CONJ</b>                                                                         | - Coniugato<br>- Conjugué<br>- Konjugat<br>- Conjugado                                                                                                                                                      | - Conjugate<br>- Conjugado<br>- Σύζευγμα                                                                                               |
| <b>MP</b>                                                                           | - Micropiastro rivestita<br>- Microplaque sensibilisée<br>- Beschichtete Mikrotiterplatte<br>- Microplaca revestida                                                                                         | - Coated microtiter plate<br>- Microplaca sensibilizada<br>- Επικαλυμμένη μικροπλάκα                                                   |
| <b>WASHB 50x</b>                                                                    | - Tampone di lavaggio<br>- Tampon de Lavage<br>- Waschpuffer<br>- Solução de lavagem                                                                                                                        | - Wash buffer<br>- Solución de lavado<br>- Ρυθμιστικό διάλυμα πλύσης                                                                   |
| <b>SUB</b>                                                                          | - Tampone substrato<br>- Substrat<br>- Substratpuffer<br>- Substrato                                                                                                                                        | - Substrate buffer<br>- Tampón sustrato<br>- Ρυθμιστικό διάλυμα υποστρώματος                                                           |
| <b>STOP</b>                                                                         | - Reagente bloccante<br>- Solution d'Arrêt<br>- Stopreagenz<br>- Solução de paragem                                                                                                                         | - Stop solution<br>- Solución de parada<br>- Αντιδραστήριο διακοπής αντίδρασης                                                         |
| <b>SB 5x</b>                                                                        | - Tampone campione<br>- Tampon Echantillons<br>- Probenpuffer<br>- Diluente de amostra                                                                                                                      | - Sample buffer<br>- Tampón Muestras<br>- Ρυθμιστικό διάλυμα δειγμάτων                                                                 |