INTENDED USE

AESKUCARE® Food Intolerance is an in vitro semi-quantitative enzyme immunoassay for the simultaneous measurement of food specific IgG₄ antibodies to 40 food extracts/extract mixtures in heparinized or Na-EDTA treated human capillary whole blood, venous whole blood, plasma or serum. It is intended for use as a point of care tool to support the diagnosis of patients presenting with various clinical symptoms associated with food intolerances. It is to be carried out by healthcare professionals experienced in the use of in vitro diagnostic methods.

PRINCIPLE OF THE PROCEDURE

AESKUCARE® Food Intolerance is a modified ELISA (Enzyme Linked Immunosorbent Assay) system. The device contains 4 rows with 5 test areas (windows) each. In each test area one standard and two specific food extracts or food mixtures are bound to the membrane. Specific IgG₄ antibodies from within the sample bind to the corresponding antigens. The bound human IgG₄ is detected by specific anti-human-IgG₄ antibodies. A secondary antibody coupled to alkaline phosphatase is added for signal amplification. After adding the color reagent bound IgG₄ is identified by an indigo colored precipitate. Food specific IgG₄ (slgG₄) levels can be read out semi-quantitatively by the naked eye for all 40 food extracts/food mixtures. The slgG₄ levels are classified into 3 levels. Read-out is performed by comparing the color saturation (intensity) of the signal lines with the one of the corresponding standard in each test field.

Precautions: 1. For human in vitro diagnostic use only. 2. Please read the entire contents of these Instructions for Use before performing the test. 3. Do not use reagents beyond their expiration dates. 4. In case of damage to the packaging please inspect the protecting bag of the device carefully for damage. Make sure the reagents are not damaged or open. In case of doubt do not use the test kit to avoid incorrect results or wrong diagnosis. 5. It is not recommended to pool any reagents. 6. The color reagent contains 5-bromo-4-chloro-3'-indolylphosphate and nitro-blue tetrazolium chloride. Avoid contact with skin. Wear suitable gloves. For more information see Material Safety Data Sheet (available on request). 7. All components are single use only. 8. Test results will not be valid if these precautions are not followed. 9. Fluid paths can be blocked by coagulated blood.

STORAGE

Store in a dark and cool place at 2-8°C/36-46°F. All reagents are ready-to-use, single dose and packaged in screw tubes.

EXPIRATION

The expiration date is printed on labels placed on the test packaging box. The expiration date of the kit is valid for all kit components, even if expiry of single components is different! After expiry all test components have to be discarded.

REAGENTS AND MATERIALS

1. **DILUT - Sample dilutor**, 350µl, 1x (orange lid) 2. **WASH - Washing solution**, 1ml each, 3x (blue lid) 3. **Ab1G4 - Antibody solution 1G4**, 800µl, 1x (grey lid) 4. **Ab2 - Antibody solution 2**, 800µl, 1x (yellow lid) 5. **COLR - Color reagent**, 800µl each (contains 5-bromo-4-chloro-3'-indolylphosphate/nitro-blue tetrazolium chloride), 2x (white lid) 6. **STOP - Stop buffer**, 1ml, 1x (green lid). **Note:** All reagents are ready for use; store at 2-8°C until expiration date.

Human material used in the AESKUCARE® Food Intolerance device components was tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.

Materials provided: instructions for use, evaluation sheet, promo card, AESKUCARE® Food Intolerance assay device, reagent rack containing 9 color coded vials, blood sampling pouch*
* **CE** certified material: disinfection swab, blood sampling capillary 50 µl with lithium heparin, swab, band-aid
* lancet: sterile, **CE** certified product, which is used as a lancing device for blood collection
* syringe (fluid application)

Additional materials required but not included in the kit: Disposable gloves, Timer

For venous whole blood and/or plasma samples: Normal venipuncture supplies, Heparinized blood collection tubes, Laboratory pipette

List of food extracts: Refer to result evaluation leaflet included in the box


TEST SAMPLES: Capillary blood from a finger prick, from earlobes or heels (in case of infants) can be used in this test. The blood sampling capillary provided in the kit should be used to collect the correct volume (50 µl). Heparinized (sodium or lithium heparin) venous whole blood, plasma, or serum can also be used as specimen in this test. The volume to be added to the sample dilutor should be 50 µl for whole blood and 20 µl for plasma or serum. In case of plasma or serum, laboratory pipettes should be used to dispense the required volume. Once added to the sample dilutor, all samples may be stored at 2-8°C for up to 72 hours. If you want to store the sample please make sure to label the test kit with the patient ID.

Note! Do not store the blood within the capillaries.

TEST PROCEDURE

Preparation: All reagents and the assay device must be adjusted to room temperature (18-30°C) for at least 30 minutes prior to use 1. Open the kit and take out the device, the syringe, the reagent rack and the blood sampling pouch. Place them on a flat surface 2. Open and assemble the reagent rack 3. Label the device with the patient ID in the grey label area.

Collection of capillary blood sample: Open the blood sampling pouch and take out its contents: alcohol wipe, safety lancet, blood sampling 50 µl capillary, gauze and band-aid. Place all items on a clean, flat surface. 1. Pinch fingertip or put it into warm water for a few moments before puncturing so as to optimize blood flow. 2. Clean the site to be pricked (ideally finger 3 or 4, lateral) with the alcohol wipe. 3. Unlock the safety lancet as shown (step A) and pull off the protective cap. 4. Place the safety lancet on the finger side to be pricked (step B) and press the button firmly.

Note! Make sure the hand is below the heart level. If the blood flow is initially slow, gently massage the finger and check if the hand is below (!) the level of heart.

Touch the drops of blood with the tip of the blood sampling capillary. Hold the blood sampling capillary horizontal or at a slight downward angle and let the blood flow until it is properly filled. Make sure not to cover the (air) hole in the plunger. The blood stop sets the sample volume (50 µl). Cover the wound with the band-aid. Immediately open the sample dilutor vial (orange lid) and put the tip of the blood sampling capillary in a position just above the vial opening. Cover the hole at the end of the plunger with your finger. Gently press down the plunger to add the blood into the sample dilutor vial. Close the sample dilutor vial and mix gently by hand or vortex briefly.

Note! Do not store the blood without having added it to the sample dilutor. Once added to the sample dilutor, all samples may be stored at 2-8°C for up to 72 hours. If untreated blood is used (i.e. without adding whole blood to sample dilutor) then the tests results will not be valid.

Collection of venous blood sample, plasma sample, or serum sample: Heparinized (sodium or lithium heparin) venous whole blood, plasma, or serum can also be used as specimen in this test. 1. Use standard laboratory techniques for collection of venous (heparinized) blood samples. 2. Spin down using standard laboratory technique to separate plasma from blood cells if desired. The volume to be added to the sample dilutor vial (orange lid) should be 50 µl for whole blood and 20 µl for plasma or serum. Laboratory pipettes should be used to dispense the required volume. Mix gently by hand after adding the sample to the sample dilutor vial.

Adding the sample: 1. Gloves should be used while performing the test 2. Open the sample dilutor vial (orange lid) containing the sample and completely withdraw its content with the syringe. Make sure no large air bubbles are formed inside. Place the syringe firmly on the fluid port of the assay device and inject the entire content of the syringe by pressing down the plunger. **Incubate for 8 minutes.** **Note!** If little backflow of the blood sample mix occurs, this does not influence the test results as long as the membrane is soaked with the blood sample mix.

Primary antibody incubation: 1. Open the first wash solution vial (blue lid), withdraw the solution completely and inject as described above 2. Open the primary antibody vial (grey lid), withdraw the solution completely and inject as described above. **Incubate for 8 minutes.**

Secondary antibody incubation: 1. Open the second wash solution vial (blue lid), withdraw the solution completely and inject as described above 2. Open the secondary antibody vial (yellow lid), withdraw the solution completely and inject as described above. **Incubate for 8 minutes.**

Color reaction: 1. Open the third wash solution vial (blue lid), withdraw the solution completely and inject as described above. 2. Open the two color reagent vials (brown tubes, white lid), withdraw the solution of the first vial completely and inject as described above. Start the timer. 3. **Within 30 sec:** Withdraw the solution of the second vial completely and inject as described above. **Incubate for 5 minutes** (total incubation time of the two steps).

Stopping the reaction: Open the stop buffer vial (green lid), withdraw the solution completely and inject as described above.

TEST PROCEDURE FOR TWO OR MORE PARALLEL TESTS

Experienced users of AESKUCARE® Food Intolerance tests may perform more than one test in parallel. It is recommended to start the tests in intervals of at least 1 minute or multiples thereof to ensure proper handling during the development steps. It is not recommended to perform more than five tests in parallel.

ASSAY RESULT READING

To evaluate the results, position the assay device with fluid port (sample inlet) pointing to the left. In this orientation the numbers of food extracts/food mixtures are positioned above and under the corresponding test areas. Results must be read in good light. Be cautious with shadows caused by the front label and avoid confusion with test signals. Start from the top left and then read out each test area until the right end of the row. Repeat row by row. This corresponds to the food list layout in the evaluation sheet. To estimate the level (1-3) of a signal its intensity must be compared with the intensity of the standard on the left side of the corresponding test area. Results are recorded in the corresponding evaluation sheet.

Note! 1. It is recommended to compare the test results with the assay result examples as shown in the illustrated Instructions for Use leaflet 2. Take care to use only the evaluation sheet that is part of the corresponding AESKUCARE® Food Intolerance kit 3. The reading of the test results should be done immediately after the test procedure to avoid confusion with other tests. Test results remain stable for at least 24 h if stored in the dark. The background color of the test areas might change with time. In case small air bubbles are present they may be removed by carefully inflating the device with some air using the syringe provided until the test areas are bubble free. If the standard band is not visible, test results are invalid.

Possible assay evaluation results

Level 1	Signal intensity less than the standard	
Level 2	Signal intensity approximately the same as the standard	
Level 3	Signal intensity greater than the standard	

Test area example		Correlation between AESKUCARE® Food Intolerance levels and food specific IgG ₄ concentration	
Standard	Signal	AESKUCARE® Food Intolerance level	slgG ₄ concentration*
		1	low
		2	moderate
		3	high
* refer to the section "interpretation of results"			

Analytical specificity: AESKUCARE® Food Intolerance detects specific IgG₄ antibodies.

Cross-reactivity with other immunoglobulin species is not expected.

Diagnostic sensitivity/specificity: The test results of AESKUCARE® Food Intolerance have been compared to slgG₄ ELISA laboratory tests of the AESKU.Group.

Restrictions and interferences: None known.

QUALITY CONTROL

It is good laboratory practice to record lot numbers of all components used.

INTERNAL CONTROLS

Each test area contains a standard (left band). This shows up in the development step as a band with a medium indigo color saturation. If the standard is missing, the test results are not valid for this test area and must either be ignored or the test must be repeated using a new test kit.

INTERPRETATION OF RESULTS

In slgG₄ diagnostics it is not usual to consider particular values as strictly defined limits. The level of slgG₄ concentration is individual specific and cannot be classified into generalized limits. The standard in AESKUCARE® Food Intolerance corresponds to a medium slgG₄ level (Level 2). Signal intensities equal or above the standard (Level 2 and 3) correspond to a moderate or high slgG₄ level. "In order to diagnose a patient, the test results must always be evaluated together with the medical history and symptoms (anamnesis).

WASTE DISPOSAL

Precautions for the handling of human blood and safe waste disposal of used parts (containing blood) must be observed. Please refer to all local/national guidelines for bio hazard handling. The rest of the kit may not be considered biohazard and can be disposed accordingly.

LIMITATIONS OF THE PROCEDURE

1. A definitive clinical diagnosis of food intolerance should not be based only on the results of a single diagnostic method. In vitro evidence of slgG₄ should always be accompanied by full medical history and analysis of symptoms 2. Each patient reacts individually, thus identical results in the test do not automatically imply the same diagnosis. Different foods with similar molecular structures or epitopes may trigger weak or severe cross-reactions. Cross-reactions must always be considered 3. Negative in vitro results may occasionally occur in patients with food intolerance symptoms that clearly correlate with food contact. 4. Level 2 or 3 results indicate the presence of slgG₄ antibodies in the patient's sample for one or more of the forty food extracts/food mixtures. High slgG₄ levels against foods other than those in AESKUCARE® Food Intolerance cannot be excluded 5. Patients with clinical symptoms suggestive of food intolerance but at a level 1 result should be referred to a specialist for further investigation 6. The binding capacity for slgG₄ antibodies may vary from food to food. Therefore, identical results for different foods may not necessarily apply exactly each slgG₄ antibody levels 7. Reliable and reproducible results will be obtained when the assay is performed according to the procedural instructions and with adherence to good laboratory practices.

EXPECTED VALUES

In general the higher the level of the test result, the higher is the concentration of food specific IgG₄. As is true with all food specific IgG₄ testing there is a potential for cross-reactivity between related and unrelated foods, which may contain similar or homologous molecules (antigens). The closer the biological relationship between different species, the greater is the degree of structural and immunological similarity of the epitopes present in the food. Accordingly, a patient who is clinically reactive to one species will likely be reactive to the other, closely related species due to immunological cross-reactivity of structurally related epitopes. On the other hand, cross-reactivity may occur also between biologically distantly related species. Some protein families are widely distributed and seem to contain highly conserved structures that can serve as similar epitopes.

WARRANTY

The performance data presented here was obtained using the procedure indicated. Any change or modification to the procedure not recommended here may affect the results, in which case Aesku.Diagnostics GmbH & Co. KG disclaims all warranties expressed, implied or statutory, including the implied warranty of merchantability and fitness for use. Aesku.Diagnostics GmbH & Co. KG and its authorized distributors, in such event, shall not be liable for damages indirect or consequential.

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ADDRESS
Manufactured by Aesku.Diagnostics GmbH & Co. KG, Mikroforum Ring 2, 55234 Wendelsheim, Germany, Tel: +49 6734 - 9622-0, Fax: +49 6734 - 9622-2222, Web: www.aesku.com
Additional company and product information can be requested at any time.

* In contrast to the situation for type 1 allergies (immediate type) there is currently no uniform evaluation system for classifying food intolerances by means of slgG₄ or slgG concentration.

HISTORY: Update in section "Material provided" and illustration of blood collection.

The instructions for use apply to AESKUCARE® Food Intolerance ONE (REF: 830101-LV2)

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