

CONTROL SERA

Positive Control Serum

Vial containing 0.5 mL positive IgG or positive IgM human control serum with protein stabilizer and 0.005% thimerosal and is prepared to give a 2+ or greater fluorescence when used as supplied (undiluted) in the AESKU.BION Immunofluorescent Test Systems. IgM Controls sold when available.

Negative Control Serum

Vial containing 0.5 mL negative human control serum with protein stabilizer and 0.005% thimerosal and is prepared to give less than a 1+ fluorescence when used as supplied (undiluted) in the AESKU.BION Immunofluorescent Test Systems.

Product Code	Product Name	Description
BBN-6110	Borrelia burgdorferi	Negative Control
BBG-6120		Positive Control, IgG
BBM-6130		Positive Control, IgM
CHN-4110	Chlamydia	Negative Control
CHG-4120		Positive Control, IgG
CHM-4130		Positive Control, IgM
CPN-4210	Chlamydia pneumoniae (TWAR) (Research Use Only)	Negative Control
CPG-4220		Positive Control, IgG
CPM-4230		Positive Control, IgM
CMN-2010	Cytomegalovirus	Negative Control
CMG-2020		Positive Control, IgG
CMM-2030		Positive Control, IgM

Product Code	Product Name	Description
EBN-5010	EBV, Early Antigen-Diffuse (Research Use Only)	Negative Control
EAD-5440		Positive Control, IgG
EBN-5010	Epstein Barr Virus VCA	Negative Control
EBG-5020		Positive Control, IgG
EBM-5030		Positive Control, IgM
EBN-5010	EBV Early Antigen - Restricted (Research Use Only)	Negative Control
EAR-5340		Positive Control, IgG
EBN-5010	EBV Early Antigen - Diffuse (Research Use Only)	Negative Control
EAD-5440		Positive Control, IgG
EBN-5010	EBV Nuclear Antigen (EBNA)	Negative Control
NA-5140		Positive Control, IgG
HSN-3510	Herpes Simplex Virus Types 1 & 2	Negative Control
HSG-3520		Positive Control, IgG
HSM-3530		Positive Control, IgM
IAN-1110	Influenza A Virus	Negative Control
IAG-1120		Positive Control, IgG
IBN-2110	Influenza B Virus	Negative Control
IBG-2120		Positive Control, IgG
MEN-6010	Measles (Rubeola)	Negative Control
MEG-6020		Positive Control, IgG
MEM-6030		Positive Control, IgM
MUN-8010	Mumps Virus	Negative Control
MUG-8020		Positive Control, IgG
MUM-8030		Positive Control, IgM
MPN-1210	Mycoplasma pneumoniae	Negative Control
MPG-1220		Positive Control, IgG
MPM-1230		Positive Control, IgM
RSN-9010	Respiratory Syncytial Virus	Negative Control
RSG-9020		Positive Control, IgG
VZN-7010	Varicella-Zoster Virus	Negative Control
VZG-7020		Positive Control, IgG
VZM-7030		Positive Control, IgM